

DIRECTA

Sendoline®

SendoCem MTA

SendoCem MTA

Root canal repair material

PRODUCT DESCRIPTION

SENDOCEM MTA is a root canal filling material conforming to ISO 6876, in a pre-loaded syringe that does not require any pre-mixing and is set by absorbing moisture from the root canal environment. It is designed to prevent irritation of pulp and secondary infections.

INDICATIONS FOR USE

The product is used for vital tooth therapy (pulp capping, Pulpotomy), root end filling, and repair of perforation.

PHYSICAL, AND CHEMICAL PROPERTIES

Setting time: 3–5 min
Radiopacity: Equivalent to 3 mm of aluminum or more in thickness
Solubility: < 3%

PRECAUTIONS FOR USE

1. WARNINGS

- Do not use in patients with allergic reactions to chemicals (specifically to Dimethyl Sulfoxide).
- For dentists' use only.
- Make sure that the product does not pass over to the mandibular canal or maxillary sinus.
- Do not use large amounts to the root canal because it is difficult to remove.

2. GENERAL PRECAUTIONS

- Check the product and packaging for abnormal condition or damage before use.
- Do not use other than for intended use.
- Keep in pouch provided before use to prevent a potential adverse reaction by carbon dioxide and moisture present in the atmosphere.
- Do not contact unhardened product with the skin or mucous membrane. If contact occurs, rinse with large amounts of water immediately.
- Wear protective gloves, mask and glasses.
- The product must be stored at room temperature avoiding high temperature and humidity.

3. ADDITIONAL PRECAUTIONS

- Make sure that the product comes out smoothly through the dedicated tip.
- Wash the treatment area with distilled water before application.
- The product shall be applied with a thickness of at least 3 mm.
- The product should be used immediately after replacing the dedicated tip because it can harden when in contact with air.

DIRECTION FOR USE

1. PREPARATIONS PRIOR TO USE

- Check the product and packaging for abnormal condition or damage before use.
- The dedicated tip should be changed with a new one prior to use. Make sure that the product comes out smoothly.

2. DIRECTION FOR USE

1) Indirect pulp capping

- Isolate the tooth with rubber dam and form cavity prep using a high speed bur.
- Place a new dedicated tip on the syringe, check if the product can be extruded smoothly.
- The product is applied to the deep dentin close to the pulp chamber, and distilled water is applied.
- After 3-minutes, remove moisture with a sterile cotton pellet, and gently pack it to avoid dead space.
- Make temporary restoration with self-curing glass ionomer etc.
- Expose radiographs for verification of achievements.

2) Direct pulp capping

- Isolate the tooth with rubber dam and form cavity prep using a high speed bur.
- When the pulp horns is exposed, it is sufficiently disinfected with undiluted NaOCl.
- Wash out with distilled water.
- Place a new dedicated tip on the syringe, check if the product can be extruded smoothly.
- Apply the product to the exposed pulp horns and control bleeding by pressing it with a cotton pellet.
- Wash out the blood with distilled water.
- If necessary, the product is applied additionally with a thickness of 3mm or more.
- Apply distilled water on the product and wait for 3-minutes.
- Remove moisture with a sterile cotton pellet, and gently pack it to avoid dead space.
- Make temporary restoration with self-curing glass ionomer etc.
- Expose radiographs for verification of achievements.

3) Pulpotomy

- Isolate the tooth with rubber dam and form cavity prep using a high speed bur.
- Remove pulp horns to the appropriate depth.
- Disinfect the surface of pulp horns with undiluted NaOCl, etc.
- Rinse the pulp horns and cavity disinfected with distilled water.
- Place a new dedicated tip on the syringe, check if the product can be extruded smoothly.
- The product is applied to the pulp horns and the surrounding hard tissue.
- Apply distilled water on the product and gently pack it with a wet cotton.
- Wash out the blood with distilled water.
- The product is applied additionally with a thickness of 3mm.
- Apply distilled water on the product and wait for 3-minutes.
- Pack the product, dry it well, and make temporary restoration with glass ionomer etc.
- Expose radiographs for verification of achievements.

4) Repair of perforation

- Isolate the tooth with rubber dam and disinfect the area.
- Dry the treatment area.
- Place a new dedicated tip on the syringe, check if the product can be extruded smoothly.
- Apply the product to the treatment area.
- Expose radiographs for verification of achievements. If the perforated site does not get filled completely, pack with an appropriate instrument, and after adding the product, a radiograph is taken.
- Perform endodontic treatment.

5) Root-end filling

- After approaching the apical area, resect apical 3 mm with surgical bur.
- Form a 3mm deep cavity on the resected apex.
- The bleeding area is compressed with a cotton pellet, etc. along with hemostatic agent.
- Discolor with methylene blue to check for cracks or accessory canals in the apical region and observe with microscope.
- Place a new dedicated tip on the syringe, check if the product can be extruded smoothly.
- Wash the reverse filling cavity with distilled water and apply the product.
- After applying the wet cotton to the product, wait 3 minutes.
- Pack the product with and appropriate instrument. If necessary, apply additional products.
- Expose radiographs for verification of achievements.

INSTRUCTIONS FOR STORAGE

- Store at room temperature avoiding high temperature and humidity.
- The new dedicated tip must be used in next use.
- After use do not recap the used dedicated tip, twist the storage cap provided with the product tightly onto the syringe head.
- Once the gray cap, which was initially fastened to the syringe, is separated, it will not be sealed even if it is fastened again. It must be sealed using the enclosed storage cap.
- Place syringe with the storage cap in the provided pouch. Seal the pouch completely.
- Tighten the storage cap immediately after use to minimize the product from contamination and infection during storage.

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Directa AB
KTR Europe GmbH