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## Manufacturer Information on Reprocessing of Medical Devices according to DIN EN ISO 17664:2004

 Manufacturer: HÄGER & WERKEN GmbH & Co. KG, Ackerstraße 1, 47269 Duisburg, Germany  
T +49 (203) 99269-0, F +49 (203) 299283, www.hagerwerken.de

Product: Fit Kofferdam Klammern

Updated: 08/2013

 355 320 – 355 321    355 338    355 340 – 355 405

 355 722 – 355 729    355 750 – 355 759    355 770 – 355 786

**Warning notices:** No special notices

**Limitation of reprocessing:** Frequent reprocessing has only got a small impact on this instrument. The instrument's shelf-life is mainly determined by wear and damage by use.

**Place of handling:** Clean surface with a disposable cloth or paper towel.

**Storage and transport:** No special requirements.

**Preparations for cleaning:** N/A

**Cleaning:** Only use cleansers according to the manufacturer's instructions. After cleaning, remove cleanser thoroughly under running water.

**Disinfection:** Use disinfecting solution according to the manufacturer's instructions. Remove disinfecting solution thoroughly under running water. Do not exceed 93 °C for mechanical disinfection.

**Maintenance:** No special requirements.

**Packaging:** Standardised packaging material may and should be used for sterilisation processes.

**Sterilisation:** Vapour sterilisation at 134 °C for 5 minutes, vapour sterilisation at 121 °C for 15 minutes

**Control/Functional check:** Sight check on damages, wear, deformation.

**Storage:** No special requirements.

The a.m. instructions have been validated as SUITABLE for preparation of a medical device and its re-use by the medical device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing site. Normally, validation and routine controls are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

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Fabricante: HAGER & WERKEN GmbH & Co. KG, Ackerstraße 1, 47269 Duisburg, Germany