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Tab. 3

## Table of Contents

|          |                                                               |           |           |                                              |           |
|----------|---------------------------------------------------------------|-----------|-----------|----------------------------------------------|-----------|
| <b>1</b> | <b>Explanation of Symbols</b> .....                           | <b>35</b> | <b>8</b>  | <b>Cleaning and Maintenance</b> .....        | <b>58</b> |
|          |                                                               |           | 8.1       | Milling machine .....                        | 58        |
|          |                                                               |           | 8.2       | Supply unit .....                            | 62        |
| <b>2</b> | <b>General Safety Instructions</b> .....                      | <b>36</b> |           |                                              |           |
| <b>3</b> | <b>Suitable Personnel</b> .....                               | <b>36</b> | <b>9</b>  | <b>Malfunctions, Repairs and Warranty</b> .. | <b>62</b> |
|          |                                                               |           | 9.1       | Malfunctions .....                           | 62        |
|          |                                                               |           | 9.2       | Repairs .....                                | 62        |
|          |                                                               |           | 9.3       | Warranty .....                               | 62        |
| <b>4</b> | <b>Machine Specifications</b> .....                           | <b>37</b> | <b>10</b> | <b>Environmental Protection</b> .....        | <b>62</b> |
| 4.1      | Milling machine .....                                         | 37        |           |                                              |           |
| 4.2      | Accessories .....                                             | 41        | <b>11</b> | <b>Technical Data and Accessories/Spare</b>  | <b>63</b> |
|          |                                                               |           |           | <b>Parts</b> .....                           | <b>63</b> |
| <b>5</b> | <b>Installation</b> .....                                     | <b>43</b> |           |                                              |           |
| 5.1      | Milling machine .....                                         | 43        |           |                                              |           |
| 5.2      | Setting -up the supply unit .....                             | 44        |           |                                              |           |
| 5.3      | Maintenance unit .....                                        | 45        |           |                                              |           |
| 5.4      | Installation of the CAM and milling<br>machine software ..... | 46        |           |                                              |           |
| 5.5      | Establishing the Network<br>Connection .....                  | 47        |           |                                              |           |
| 5.6      | Starting-up the milling machine ..                            | 48        |           |                                              |           |
| <b>6</b> | <b>Application and Operation</b> .....                        | <b>49</b> |           |                                              |           |
| 6.1      | Switching between wet and dry<br>processing .....             | 50        |           |                                              |           |
| 6.2      | Changing the blank holder .....                               | 52        |           |                                              |           |
| 6.3      | Program execution .....                                       | 53        |           |                                              |           |
| 6.4      | Checking the workpiece during a<br>running program .....      | 53        |           |                                              |           |
| <b>7</b> | <b>Handling the Coolant/Lubricant (C/L)</b> .                 | <b>54</b> |           |                                              |           |
| 7.1      | C/L analysis .....                                            | 54        |           |                                              |           |
| 7.2      | Removing the supply unit from the<br>set-up cart .....        | 55        |           |                                              |           |
| 7.3      | Emptying the C/L container .....                              | 56        |           |                                              |           |
| 7.4      | Cleaning the C/L container .....                              | 56        |           |                                              |           |
| 7.5      | Filling in C/L .....                                          | 57        |           |                                              |           |
| 7.6      | Installing the supply unit into the<br>set-up cart .....      | 57        |           |                                              |           |
| 7.7      | Cleaning the C/L system .....                                 | 57        |           |                                              |           |

# 1 Explanation of Symbols

## Warning indications

Warning indications in the text are marked with a colour-backed triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- \_ **NOTE** means that property damage can occur.
- \_ **CAUTION** means that light to fairly serious personal injury can occur.
- \_ **WARNING** means that serious personal injury can occur.
- \_ **DANGER** means that serious personal injury can occur.

## Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

## Other symbols in the Manual

| Symbol | Meaning                                                                  |
|--------|--------------------------------------------------------------------------|
| >      | Item of an operation description                                         |
| _      | Item of a list                                                           |
| •      | Subitem of an operation description or a list                            |
| [3]    | Numbers in square brackets refer to position numbers in graphics/figures |

## Other symbols on the machine

| Symbol | Meaning                               |
|--------|---------------------------------------|
|        | USB port (socket type B)              |
|        | Network connection (Ethernet)         |
|        | Slot for SD card                      |
|        | Control cable connection, extraction  |
|        | Fuse                                  |
|        | Compressed-air connection             |
|        | Control cable connection, supply unit |
|        | Supply connection, coolant/lubricant  |

## 2 General Safety Instructions

When installing, starting-up and operating the machine, always observe the following safety instructions:

### CAUTION:

A defective machine can lead to malfunctions! On detection of damage or a functional defect of the machine:

- ▷ Label the machine as defective.
- ▷ Prevent further operation until the machine has been repaired.

### NOTE:

Possible machine damage through insufficient extraction during dry processing!

- ▷ Operate the machine only with the extraction system recommended by the manufacturer or another compatible extraction system (see page 46).

### NOTE:

Possible machine damage through insufficient cooling during grinding operations!

- ▷ For grinding operations, operate the machine only with the supply unit specifically intended for it.

### NOTE:

Damage through escaping coolant/lubricant!

- ▷ Regularly check the coolant/lubricant connection for tightness against leaks.

### NOTE:

- ▷ Switch the machine off when not in use or unsupervised for longer periods, e.g., overnight. This method of saving electrical energy also benefits the environment.

## 3 Suitable Personnel

### NOTE:

Starting-up and operation of the machine may only be carried out by trained specialised personnel.

## 4 Machine Specifications

### 4.1 Milling machine

This User Manual describes different models of the milling machine. Therefore, the representations can deviate from the actual appearance of the machine.

#### 4.1.1 Delivery Scope

- \_ Digital milling machine
- \_ Mains cable
- \_ USB cable (5m)
- \_ Network cable (5m)
- \_ Control cable, extraction
- \_ Service unit including compressed-air hoses
- \_ Service set for spindle
- \_ Torque wrench for collet
- \_ CAM software
- \_ Auto-calibration set
  - Calibration blank
  - Calibration pin
- \_ Set for manual calibration
  - Test blank
  - Cutter 2.5 mm
  - Slide caliper
- \_ Torque wrench (4 Nm) for blank holder
- \_ Torque wrench (3 Nm) for fastening blank
- \_ Blank holder 98 mm
- \_ Extractor adapter
- \_ Extraction hood
- \_ Drain sieve
- \_ Cleaning brush
- \_ Cleaning brush

The milling machine is supplied with an automatic calibration set.

To check the zero point, a test milling run was carried out at the factory prior to delivery. Milling chips may have remained from this. These are no reason for rejection, as the check process contributes significantly to the accuracy of the machine.

- ▷ After unpacking, check the machine for completeness and possible transport damages. Please claim any transport damages immediately with your supplier.

### 4.1.2 Intended Use

The milling machine is a PC-controlled milling and grinding machine for the fabrication of dentures by means of:

- \_ Dry processing of blanks made of
  - pre-sintered zirconium oxide
  - unsintered non-precious metals (NPMs)
  - Wax
  - Plastics
- \_ Wet processing of blanks made of
  - Sintered glass ceramics
  - Plastics

Using blanks and tools not approved by the manufacturer can damage the machine and make the product unusable. For such cases, the manufacturer shall assume no liability whatsoever.

Unauthorised modifications/alterations of the machine shall also void the warranty.

A compressed-air connection is required for operation of milling machine (see page 45).

Dry processing is permitted only in conjunction with an appropriate extraction system (for technical requirements, see page 46).

Wet milling and grinding is permitted only in conjunction with the supply unit (accessory).

### 4.1.3 CE Declaration of Conformity

This product was designed and manufactured based on careful selection of the harmonised standards to be observed, as well as additional technical specifications. It thus corresponds with the state-of-the-art and ensures maximum safety.

In terms of design and performance, this product complies with the European and UK Directives and the supplementary national requirements. Conformity has been confirmed with the CE and UKCA marking.

Directives:

- \_ 2011/65/EU:2011-06-08 RoHS II Directive
- \_ 2006/42/EG:2006-05-17 Machinery Directive
- \_ 2014/30/EU:2014-02-26 EMC Directive
- \_ S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016
- \_ S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

Applied harmonised standards:

- \_ DIN EN ISO 12100:2011-03
- \_ DIN EN 63000:2019-05
- \_ DIN EN 61010-1:2020-03
- \_ DIN EN 55011:2022-05
- \_ DIN EN 61000-4-4:2013-04
- \_ DIN EN 61000-4-5:2019-03
- \_ DIN EN ISO 16090-1:2019-12

The product's declaration of conformity is available upon request from the manufacturer.

## 4.1.4 Components and Interfaces


Fig. 1 Machine overview, front

- |                                 |                        |
|---------------------------------|------------------------|
| 1 Touch-display                 | 6 C/L drain            |
| 2 Measuring device, tool length | 7 Tool holder          |
| 3 C/L nozzles                   | 8 RFID read/write unit |
| 4 Milling spindle               | 9 Blank holder         |
| 5 Extraction hood               |                        |



Fig. 2 Machine overview, rear

- 1 Slot for SD memory card 
- 2 Ethernet interface 
- 3 USB data port 
- 4 Control interface for extraction 
- 5 WLAN antenna
- 6 Compressed-air inlet
- 7 Compressed-air hose, 250 mm
- 8 Maintenance unit
- 9 Control interface for supply unit 
- 10 Connection, C/L supply line 
- 11 Compressed-air connection, milling spindle
- 12 Connection, C/L drain
- 13 Power supply connection
- 14 Fuse compartment
- 15 Main switch
- 16 Extraction-hose connection

## 4.2 Accessories

### 4.2.1 Milling and grinding tools

Various cutting and grinding tools are available for the milling machine (see catalogue).

### 4.2.2 Supply unit

#### Intended Use

The supply unit enables wet processing in the milling machine. It ensures the application-efficient supply, drainage and filtering of the required coolant/lubricant.

Operation with other machines than the milling machine is not permitted.

#### Delivery Scope

- \_ Supply unit with
  - Container for coolant/lubricant (C/L)
  - Filter for (C/L)
- \_ Control cable
- \_ C/L supply line
- \_ Piping kit
- \_ Swarf sieve
- \_ Splash guard
- \_ 2 litres of coolant/lubricant
- \_ Analytics set
  - Refractometer
  - Test for pH value
  - Test for nitrite content
  - Test record

## Components and Interfaces


Fig. 3 Connections of the supply unit

- 1** C/L drain cock
- 2** C/L container
- 3** Connection for C/L supply line
- 4** Control interface for supply unit
- 5** C/L drain

## 4.2.3 Set-up cart

The set-up cart is used as a mobile carriage for the milling machine. The feet of the milling machine are positioned into the retainers of the set-up cart. The supply unit can be placed into the bot-

tom compartment of the set-up cart. An extraction system can also be accommodated there, if required.


Fig. 4 Supply unit in set-up cart

- 1 Set-up cart
- 2 C/L drain
- 3 C/L supply line
- 4 Supply unit
- 5 C/L drain cock

## 4.2.4 Extraction

The extraction recommended by the manufacturer is required for dry processing with the milling machine. It is equipped with a special filter bag and a HEPA micro filter, which filters up to 99.97 % of fine dust particles (i.e. zirconium oxide dust particles) and is classified in the fine dust particle category M.

The separate socket for the control cable allows for automatic operation.

## 5 Installation

### 5.1 Milling machine

#### Countries with 100 - 115 V mains voltage

The machine is preset for a mains voltage of 230 V. When the mains voltage at the set-up location is 100 - 115 V:

- ▷ Pull out drawer ([14] in Fig. 2 on page 40) above the mains plug.
- ▷ Replace both T3.15 A fuses with the provided T6.3 A fuses.

#### 5.1.1 Setting up the milling machine

Please refer to the separate milling machine set-up instructions.

- \_ The machine is intended exclusively for use within dry, closed rooms.
- \_ A clearance of at least 200 mm is to be observed sideways, to the rear and upwards. Minimum space requirements for the milling machine (incl. connections):
  - Width: 930 mm
  - Depth: 795 mm
  - Height: 980 mm
- \_ The empty weight of the machine is between 78 kg and 82 kg. The set-up surface must withstand an appropriate load.
- \_ During operation, the room temperature must be between 18 °C and 30 °C; avoid large variations in temperature.
- ▷ Set up the machine on a sturdy work table or work bench (off the ground and clear of walls), or on the set-up cart available from the manufacturer.
- ▷ Adjust the height of the machine by screwing the feet in or out, ensuring that the machine is positioned firmly on all four feet.
- ▷ Mount the service unit (see page 45).
- ▷ Insert the extraction hose of the extraction system to the connection on the rear side of the milling machine.
- ▷ Connect the milling machine according to the connection diagram in the set-up instructions.

## Setting up the milling machine on the set-up cart

A clearance of at least 200 mm is to be observed sideways, to the rear and upwards. Space requirements for the milling machine on the set-up cart incl. connections:

- \_ Height: 1660 mm
- \_ Width: 930 mm
- \_ Depth: 860 mm

The set-up surface must be level.


Fig. 5

- 1 Milling machine
- 2 Set-up cart

### NOTE:

Possible damage of the electric lines and the C/L supply line:

- ▷ When moving the set-up cart, pay attention that no lines are driven over.

- ▷ Position the milling machine onto the set-up cart, ensuring that the legs engage into the retainers.
- ▷ Move the set-up cart with the milling machine to the set-up location.
- ▷ Lock the casters of the set-up cart with the brakes.

## 5.2 Setting up the supply unit

Set up the supply unit below the work table/bench of the milling machine or into the set-up cart available from the manufacturer.

In order for the C/L drain pipe to be mounted, the table must be set up with a certain clearance to the wall (see set-up instructions).

- ▷ Fill coolant/lubricant into the C/L container (see section 7.5 on page 57).

### NOTE:

Possible damage of the electric lines and the C/L supply line:

- ▷ When moving the set-up cart, pay attention that no lines are driven over.
- ▷ Connect the C/L supply line to the milling machine and the supply unit.
- ▷ Plug the C/L drain pipe together and connect it to the milling machine and the supply unit (see set-up instructions).
- ▷ Connect the control cable to the control interface for the supply unit on the rear side of the milling machine ([9] in Fig. 2) and to the supply unit ([4] in Fig. 3).
- ▷ Set up the C/L container or insert it into the set-up cart, as required (please refer to section 7.6).

### 5.3 Maintenance unit

The spindle is equipped with a sealing-air feature. This sealing air prevents shavings/chips and dust from entering the spindle. The service unit filters possible contamination in the sealing air and thus prevents costly machine damage.

#### Installing the service unit


Fig. 6 Maintenance unit

- 1 Compressed-air inlet
- 2 Water separator
- 3 Compressed-air hose, 250 mm
- 4 Control interface for supply unit
- 5 Connection, C/L supply line
- 6 Compressed-air connection, milling spindle
- 7 Connection, C/L drain pipe

- ▷ Mount the service unit on the rear side of the machine using the screws.
- ▷ Connect the 250 mm compressed-air hose [3] to the right-hand connection of the service unit and to the compressed-air connection [5] of the milling machine.
- ▷ Connect the 2000 mm compressed-air hose to the compressed-air inlet [1] of the service unit and the other end to the compressed-air supply of the laboratory.

#### Requirements for the compressed-air supply

For operation, the compressed air must fulfil the following minimum requirements:

- \_ Flow rate: 50 l/min (1.77 CFM)
- \_ Pressure: 6 bar / 87 psi

Requirements on air purity:

- \_ Solid contaminants: Class 3; better than 5 µm for solids
- \_ Water content: Class 4; Max. pressure dew point +3 °C
- \_ Total oil content: Class 2; Max. oil content 0.1 mg/m<sup>3</sup>

In continuous operation, the pressure and the air-flow rate are monitored by the milling machine. When the required values are fallen below, the milling program stops. Once the correct air values are available again, the milling program can be continued after fault acknowledgement in the software.

### Requirements on air purity for the spindle

When another service unit than the one supplied is to be used, the following connection conditions according to

DIN-ISO 8573-1 must be ensured:

- \_ Solid contaminants: Class 3  
(particle size max. 5 µm, particle density max. 5 mg/m<sup>3</sup>)
- \_ Water content: Class 4  
(max. pressure dew point +3 °C, water content max. 6000 mg/m<sup>3</sup>)
- \_ Total oil content: Class 2;  
(max. oil content 0.1 mg/m<sup>3</sup>)

### Recommended extraction system requirements

When another extraction system is to be used, it must meet the following requirements:

- \_ Extraction capacity: 240 m<sup>3</sup>/h
- \_ Suitability for zirconium-oxide dusts
- \_ HEPA microfilter (97.97 %), filter class H12, dust category M

The connection of the extraction hose to the machine is made with an adapter (in delivery scope). The adapter fits extraction hoses with a 38 mm interior diameter.

### 5.4 Installation of the CAM and milling machine software

- ▷ Insert the CAM software CD-ROM into the drive of the PC.  
The software to be installed is on the CD.
- ▷ Execute the “.exe” file and follow the directions on the monitor.

### Connecting the interfaces


Fig. 7 Interfaces on the rear side of the milling machine

- 1 Slot for SD memory card
- 2 Ethernet interface
- 3 USB data port
- 4 Control interface for extraction
- 5 Main switch
- 6 Power supply connection

- ▷ Make sure that the main switch [6] of the milling machine is in the 0 position.
- ▷ Check if the SD memory card is properly seated in the slot [1] on the rear side of the milling machine.
- ▷ Connect the mains cable to power supply connection [7] of the milling machine as well as to a mains socket outlet.

- ▷ Connect the control cable to the control interface for extraction [5] on the rear side of the milling machine and to the interface on the rear side of the extraction system.

For installation with supply unit:


Fig. 8 Connections of the supply unit

- 1** Control interface for supply unit
- ▷ Connect the control cable to the control interface for the supply unit on the rear side of the milling machine ([4] in Fig. 7) and to the supply unit ([1] in Fig. 8).

### Connecting to the PC

- ▷ Connect a USB cable to the USB interface ([3] in Fig. 7) of the milling machine and to a free USB port of the PC.  
The maximum permitted length of the USB cable is 5 m.
- ▷ Insert the CAM Dongle in a free USB port on the rear side of the PC.


Fig. 9 CAM Dongle

## 5.5 Establishing the Network Connection

To establish the network connection, the milling machine and the PC must first be connected via USB.

- ▷ Connect the Ethernet interface ([2] in Fig. 7) of the milling machine with the PC or with the network using a network cable.  
The maximum permitted length of the network cable to the PC or to the network router is 100 m.
- ▷ Switch the milling machine on.
- ▷ Start the milling machine software.
- ▷ After initializing the milling machine in the software under *Settings > Communication*, set the connection type from USB to Ethernet and adapt the IP address, if required (see description of the milling machine software).
- ▷ Switch the milling machine off and then on again.
- ▷ Initialize the milling machine in the software.
- ▷ After the software has detected the Ethernet connection, remove the USB cable.

When operating via a network, up to eight machines can be connected and controlled via the CAM software.  
The PC and the milling machine must both be in the same sub-network.

## 5.6 Starting-up the milling machine

Calibration of the milling machine is advisable after installation for commissioning, but not mandatory. The milling machine can adjust itself via auto-calibration. This process can be started on the display under *Maintenance*. All necessary steps are explained on the display. In the following, the calibration is described via the milling process. This is to be performed in exceptional cases and only upon support request.

The cutters can be inserted in all slots of the milling machine. By means of the integrated RFID chip, the machine automatically recognizes the cutter. When a standard tool has reached its wear limit, the milling machine automatically accesses the respective sister tool.


Fig. 10 Tool holder

- ▷ For calibration, insert the "760660 Roto RFID 2.5 Zl" cutter from the delivery scope.

- ▷ Insert the selected blank into the holder in the milling machine.


Fig. 11 Inserting a blank

- ▷ Fasten the blank with the bayonet lock.
- ▷ Lock the screws with the provided torque wrench.


Fig. 12 Tightening a blank

- ▷ Mount the respective hood.
  - Extraction hood for dry processing
  - Splash guard for wet processing

**For wet processing**

- ▷ Insert the swarf sieve with the wide side facing toward the rear.

When the chip sieve is not inserted for wet machining, chips are flushed directly into the tank. This leads to more frequent changing of the filter bag.

The swarf sieve is included in the delivery scope of the supply unit. It can remain inserted in the machine during dry processing.

- ▷ Shut the door of the milling machine.
- ▷ Switch the milling machine on via the main switch.
- ▷ Switch the PC on.
- ▷ Start the milling machine software.

When the software is started, a machine initialization takes place.

If the milling machine is not detected:

- ▷ Restart the PC.

After re-positioning, it is recommended to check the system axes. This can take place through:

- ▷ Activation of auto-calibration in the service menu of the milling machine.
- ▷ A manual test milling run (in exceptional cases and only upon support request).

**6 Application and Operation****WARNING:**

Danger from shavings/chips or tool fragments being thrown from the machine!

- ▷ Make sure that the door of the milling machine is always closed during processing!

**NOTE:**

Insufficient milling results and increased wear when dry processing!

- ▷ Do not operate the milling machine without extraction!

**NOTE:**

Insufficient milling/grinding results and increased wear when wet processing!

- ▷ Do not operate the milling machine without supply unit!

## 6.1 Switching between wet and dry processing

The milling machine can be used both for wet and dry processing.

In delivery condition, it is set up for dry processing.

### 6.1.1 Switching from dry to wet processing


Fig. 13

- 1 Extraction hood
  - 2 Lever of the extraction hood
  - 3 Swarf sieve
  - 4 Splash guard
  - 5 Lever of the splash guard
- ▷ Disengage lever [2] and pull off the white extraction hood [1] toward the front.
  - ▷ Connect the extraction hose to the extraction opening in the milling chamber.
  - ▷ Activate the extraction system via the display under *Service > Exhaustion*.
  - ▷ Vacuum out the milling chamber.

When the milling chamber is clean:

- ▷ Deactivate the extraction system via the display under *Service > Exhaustion*.
  - ▷ Pull off the extraction hose.
  - ▷ Insert the cleaned swarf sieve [3]. The sieve serves as the first coolant/lubricant stage for coarse particles.
- The wide side faces toward the rear.

For wet processing, it is mandatory to use the swarf sieve [3]. The sieve serves as the first coolant/lubricant stage for coarse particles.

- ▷ Check if the outlet holes of the C/L nozzles are unobstructed.

- ▷ Check if the door gasket is clean and free of swarf/chips.
- ▷ Check if the C/L drain pipe is free of obstructions and if the piping on the outside is securely connected.

- ▷ Attach the splash guard [4] to the extraction opening and lock it with lever [5].
- ▷ Transfer the milling file for wet processing.

## 6.1.2 Switching from wet to dry processing


Fig. 14

- 1 Extraction hood
- 2 Lever of the extraction hood
- 3 Swarf sieve
- 4 Splash guard
- 5 Lever of the splash guard

- ▷ Release lever [5] and pull off the splash guard [4] toward the front.
- ▷ Remove the swarf sieve [3].
- ▷ Clean and dry the milling chamber and the swarf sieve (see page 58).

 The swarf sieve [3] is not required for dry processing.

- ▷ Insert the blank and insert or change the tooling, as required.
- ▷ Mount the white extraction hood [1] to the extraction opening and lock it in place with lever [2].
- ▷ Check the air extraction.
- ▷ Transfer the milling file for dry processing.

## 6.2 Changing the blank holder

In delivery condition, the milling machine is equipped with a 98 blank holder. For processing varying materials, it may be required to use a different blank holder.

To change the blank holder:

- ▷ With the door closed, select *Service* > *Change blankholder* on the display of the milling machine.

The current tool is placed down. The blank holder [2] moves to the change position.


Fig. 15 Changing-position of the blank holder

**1** Fastening screw for blank holder

**2** Blank holder

- ▷ Loosen fastening screw [1].
- ▷ Depending on the model, pull off the blank holder toward the right or toward the bottom.
- ▷ Clean the locating bore and surfaces.
- ▷ Check the pin of the blank holder and clean if required.
- ▷ Insert the blank holder with the pin to the stop in the locating bore.

### NOTE:

Inaccurate working results or damage to the blank-holder fixture possible!

- ▷ Tighten fastening screw [1] with a torque of 4 Nm.

The torque wrench included in the delivery scope is appropriately set.

- ▷ Tighten fastening screw [1] with with the supplied torque wrench.

- ▷ Close the milling chamber door.

The machine moves to the starting position.

### 6.3 Program execution

Before executing a program, it must first be loaded (see CAM software). After loading, a new job appears on the milling machine display.

- ▷ Start the new job on the display of the machine with *Start job*.

The program can be stopped:

- ▷ On the display of the milling machine, press *Pause job*.

To continue the program:

- ▷ With the door of the milling chamber closed, press *Continue job* on the display.  
The program is continued at the position where it was stopped.

On completion of the program, the machine returns to the starting position.

### 6.4 Checking the workpiece during a running program

- ▷ On the display of the milling machine, press *Pause job*.  
The program is interrupted. The spindle moves to the rear.
- ▷ Open the milling chamber door.
- ▷ Visually check the workpiece.
- ▷ Close the milling chamber door.
- ▷ On the display of the milling machine, press *Continue job*.

The tool is picked up again. The spindle travels to the position, where the program was interrupted.

## 7 Handling the Coolant/Lubricant (C/L)

### CAUTION:

Health hazards are possible through improper handling of C/L!

- ▷ Use only C/L recommended by the manufacturer!
- ▷ Observe the safety instructions and handling regulations for C/L!

### 7.1 C/L analysis

For legal and technical reasons, the coolant/lubricant must be analysed weekly.

- ▷ Observe country-specific laws and guidelines on coolant/lubricant!

Regular analysis of the C/L is required for proper function of the milling machine. The results of the analyses must be documented and archived.

- ▷ Check the following coolant/lubricant properties weekly:
  - Concentration
  - pH-value
  - Nitrite content

The analysis is performed with test strips and refractometer. The test strips are easy to use and allow the determination of the pH value and the nitrite content. The concentration of the C/L is measured with the refractometer.

- ▷ Remove the blank from the holder.
- ▷ Mount the splash guard.
- ▷ Place a clean, flat receptacle in the centre of the interior below the splash guard.
- ▷ Shut the door.

- ▷ On the display of the milling machine in the *Service menu*, activate the coolant pump for 10 seconds.
- ▷ Open the door and remove the receptacle.
- ▷ Measure the pH and nitrite value:
  - Immerse the reaction zone of the pH and nitrite measuring strip into the C/L.
  - Shake off excess fluid.
  - After one minute, compare the reaction zone with the colour scale and read off the measured value.
- ▷ Measure the concentration:
  - Remove C/L with the pipette and apply it to the measuring surface (glass area).
  - Close the lid of the refractometer.
    - The measuring surface must be completely wetted with liquid.
  - Hold the refractometer against a light source and look through the ocular (eyepiece).

The reference values and the information on measures in case of exceeding or falling below the reference values can be found on an information sheet in the delivery scope of the supply unit.

## 7.2 Removing the supply unit from the set-up cart

▷ Carefully lift the supply unit [2] out of the set-up cart.

**NOTE:**  
 Possible damage of the electric lines and the C/L supply line:  
 ▷ When moving the set-up cart, pay attention that no lines are driven over.

- ▷ Move the set-up cart out from below the table.
- ▷ Remove the left-hand side wall of the set-up cart.
- ▷ Pull the C/L drain pipe [1] and the C/L supply line [3] at the top from the supply unit.
- ▷ Pull off control cable [2].


Fig. 17

- 1 Set-up cart
- 2 Supply unit


Fig. 16 Releasing the interfaces at the supply unit

- 1 C/L drain pipe
- 2 Control cable
- 3 C/L supply line

**NOTE:**  
 Possible damage to the drain cock!  
 ▷ Move/carry the supply unit only by the handles!

## 7.3 Emptying the C/L container

### NOTE:

Possible damage to the drain cock!

- ▷ Move/carry the supply unit only by the handles!

- ▷ Place the supply unit [1] with its rear side facing frontward on a table.


Fig. 18

- 1 Supply unit
- 2 Drain cock
- 3 Disposal receptacle

- ▷ Position a disposal receptacle [3] below the drain cock [2].
- ▷ If required, mount a hose to the drain cock and conduct it into the disposal receptacle.

- ▷ Open the drain cock.

The C/L flows out of the C/L container into the disposal receptacle.

When the C/L container is empty:

- ▷ Shut the drain cock again.

Dispose of the C/L in accordance with the applicable regulations and laws. Proper disposal of the C/L must be documented.

## 7.4 Cleaning the C/L container

### NOTE:

The filter and the C/L container may only be cleaned with water.

- ▷ Do not use cleaning agents.

- ▷ Unscrew and remove the four star-knob nuts on the top side of the supply unit.
- ▷ Remove the cover with the drain pipe.
- ▷ Remove the filter bag and empty it.
- ▷ Clean the filter bag.
- ▷ Rinse and clean the C/L container.

## 7.5 Filling in C/L

The water for the coolant/lubricant must have drinking-water quality. The optimal hardness range is between 10 - 20 °dH (corresponds to 1.8 - 3.6 mmol total hardness per litre). If the hardness is far beyond this range, we recommend the addition of distilled water. For grinding glass-ceramic materials, it is generally recommended to use distilled water.

- ▷ Fill 20 liters of water and 1 liter of coolant concentrate into the C/L container and mix well.
- ▷ Insert the filter bag into the C/L container.
- ▷ Reattach the cover with the drain pipe and secure it with the four star-knob nuts.

## 7.6 Installing the supply unit into the set-up cart

- ▷ Lift the supply unit into the set-up cart.
- ▷ Connect the control cable.
- ▷ Connect the return pipe and the C/L supply line at the top to the supply unit.
- ▷ Mount the left-hand side wall of the set-up cart.

### NOTE:

Danger of damaging cables and hoses!

- ▷ When moving the set-up cart, do not pinch or drive over any cables or hoses.
- ▷ Connect the C/L drain pipe, the C/L supply line and the control cable to the milling machine.

## 7.7 Cleaning the C/L system

- ▷ Remove and clean the swarf sieve.
- ▷ Clean the milling chamber.
- ▷ Reinsert the swarf sieve.
- ▷ Mount the splash guard.
- ▷ Empty the C/L container (see page 56).
- ▷ Clean the C/L container (see page 56).
- ▷ Fill 5 liters of water and 250 ml of coolant concentrate into the C/L container and mix well.
- ▷ Flush the milling machine for 15 minutes.
- ▷ Empty the C/L container.
- ▷ Insert the filter bag.
- ▷ Fill 20 liters of water and 1 l of coolant concentrate into the C/L container and mix well.
- ▷ Fully connect the supply unit again.

## 8 Cleaning and Maintenance

The milling machine has interactive cleaning workflows in the menu area *Maintenance*.

The following applies for all cleaning runs:

- \_ When due, they are automatically indicated on the milling machine display.
- \_ They can be conveniently called up and carried out, documented with pictures, step by step on the display of the milling machine.
- \_ They are documented.

### 8.1 Milling machine

#### 8.1.1 Cleaning

The milling machine must be cleaned after each working day.

- ▷ Empty the swarf sieve frequently, especially when wet grinding plastic.  
When the swarf sieve is too full, it is possible that the coolant/lubricant will not drain off.

#### Cleaning after dry processing

**NOTE:**

Zirconia dust is abrasive!

- ▷ Clean the window in the door carefully, using only a soft cloth; do not scrub or scour.

- ▷ Open the milling chamber door.
- ▷ Replace the extraction hood with the cleaning hose with brush head.
- ▷ On the display of the milling machine, activate *Exhaustion* and *Sealing air*.

- ▷ If required, remove the swarf sieve and dispose of accumulated milling dust/shavings away from the machine.

In this manner, the filter bag in the extraction system will not fill up as quickly and can be used longer.

- ▷ Thoroughly vacuum out the milling chamber with the cleaning hose.
- ▷ Carefully clean the RFID read/write unit with a clean cloth.
- ▷ Carefully clean the tool shanks.

Residual zirconia dust from dry processing can clog the filter bag in case of subsequent wet processing.

- ▷ On the display of the milling machine, deactivate *Exhaustion* and *Sealing air*.

#### Cleaning after wet processing

**NOTE:**

Possible machine damage through dried-in coolant/lubricant and bound in milling dust.

After wet processing:

- ▷ Clean and dry the milling chamber prior to downtimes lasting several hours.

**NOTE:**

The extraction is not designed for wet vacuuming. The extraction may not be used for cleaning the interior after wet processing!

- ▷ Remove shavings/chips manually.

**NOTE:**

The milling chamber may only be cleaned with water.

- ▷ Do not use cleaning agents.

- ▷ Remove the swarf sieve and dispose of accumulated milling dust/shavings.
- ▷ Clean the milling chamber and wipe it dry using a soft cloth.

 Damage caused by not removed grinding residue are not covered by the warranty.

### 8.1.2 Regular inspections

- \_ Condition of the bellow
- \_ Condition of hoses and drain pipes
- ▷ Keep the door gasket and the stopping surface on the door clean. Ensure not to damage the sealing lips.

### 8.1.3 Weekly maintenance

 **NOTE:**  
 Danger of damaging the machine!  
 ▷ Do not use compressed air, ultrasonic or a steam jet for the cleaning!

- ▷ On the display of the milling machine, select *Maintenance > Clean machine.*
- ▷ Follow the instructions on the display.

### Spindle

The provided spindle service set and the collet torque wrench are required for maintenance of the spindle.


Fig. 19 Service set for spindle

- 1 Collet grease
- 2 Collet brush
- 3 Felt cone
- 4 Ejector pin
- 5 Torque wrench for collet

- ▷ On the display of the milling machine, select *Maintenance > Clean chuck.*

In order to clean the spindle:

- ▷ Remove any dust from the cap, the C/L pipes and the spindle using a dry, soft cloth. Do not use cleaning agents!

To remove the collet:

- ▷ Apply collet torque wrench [5] to collet.


Fig. 20

- ▷ Using the collet torque wrench, unscrew the collet out of the spindle.


Fig. 21

- ▷ Clean the collet holder using the felt cone [3]. The collet holder must be free of shavings/chips and debris.


Fig. 22

- ▷ Clean the inside of the collet with the fine collet brush [2].


Fig. 23

- ▷ Apply a light coat of the collet grease ([1] in Fig. 19 on page 59) from the service set to the outer cone of the collet. This improves the sliding properties and increases the clamping force of the collet.

### NOTE:

Only use the provided collet grease for lubrication of the collet cone. Do not use other types of grease.

- ▷ Screw the cleaned and lubricated collet in reverse order to the stop back into the spindle.

### C/L nozzles

Do not clean the C/L nozzles with the collet brush!

- ▷ Inspect the C/L nozzles weekly.
- ▷ Clean the outlet holes as required.

When the outlet holes are clogged:

- ▷ Loosen the Allen screw on the side of the spindle cap.


Fig. 24

▷ Pull off the cap.


Fig. 25

▷ Blow out the holes from the inside of the cap using a steam cleaner.

**CAUTION:**  
The cap heats up!

- ▷ Apply a coat of grease to the C/L pipes.
- ▷ Clean the seat of the cap on the spindle.
- ▷ Mount the spindle cap again with the screw at the right-hand side. Pay attention to the correct connection of the C/L pipes.

### Maintenance unit

- ▷ Visually check the service unit for:
  - Visible particles, debris or particulate matter within the acrylic glass (signs for contaminated compressed air)
  - Amber discolouration of the filter element (sign of oil in the compressed air)
  - Water in the filter (sign of water in the compressed air)

Once any of the listed impairments is given:

- ▷ Correct the problem with the compressed air.
- ▷ Replace the complete service unit.

## 8.1.4 Monthly maintenance

### Checking the zero point

It is recommended to check the zero point once every month. We recommend the auto-calibration of the milling machine for this purpose.

- ▷ On the display of the milling machine, select *Maintenance > Calibrate machine.*

### Spindle

The spindle must be in the starting position.

- ▷ Manually turn the shaft of the spindle through (at least ten times).


Fig. 26 Rotating the spindle

## 8.1.5 External maintenance

External maintenance is necessary after 1000 operating hours. The milling machine software will inform you appropriately.

## 8.2 Supply unit

### 8.2.1 Inspections and maintenance measures

Regular inspections:

- \_ Condition of hoses and drain pipes
- \_ Filling level of the C/L in the supply unit

 When the filling level in the C/L container has fallen below the minimum level, the milling program is interrupted and an error message is indicated.

- ▷ Fill the C/L container (see page 57).
- ▷ Acknowledge the error in the software.

### 8.2.2 Weekly maintenance

- ▷ Analyse the C/L (see separate analysis instructions with test protocol).

### 8.2.3 Half-yearly maintenance

- ▷ Clean the C/L system (see page 56).

 The half-yearly cleaning is indicated by the program. Upon activation of the system cleaning, the indication is blanked out again.

## 9 Malfunctions, Repairs and Warranty

### 9.1 Malfunctions

In case of malfunctions:

- ▷ Restart the software.
- ▷ Restart the milling machine.
- ▷ Restart the PC.

### 9.2 Repairs

Repairs may be carried out only by trained, specialised personnel.

### 9.3 Warranty

The warranty complies with the statutory provisions. For more information, please refer to our General Terms of Business.

## 10 Environmental Protection

### Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally friendly and recyclable.

### Waste equipment

Waste equipment contain materials that should be sorted for environmental-friendly recycling.

- ▷ After the service life has elapsed, dispose of the machine via the public disposal systems in an environmentally-friendly manner.

The plastic parts are labelled. Thus, they can be sorted for recycling or disposal.

### Coolant/Lubricant (C/L)

- ▷ Dispose of the C/L in accordance with the applicable regulations and laws.

Proper disposal of the C/L must be documented.

11 Technical Data and Accessories/Spare Parts

 Subject to changes.

**Technical data of the milling machine**

|                                         | Unit  | Value                                                                                      |
|-----------------------------------------|-------|--------------------------------------------------------------------------------------------|
| Art. No.                                | -     | 181350                                                                                     |
| Dimensions (D x W x H)                  | mm    | 595 x 530 x 780                                                                            |
| Weight                                  | kg    | 78 - 82                                                                                    |
| Electrical connections                  | V/Hz  | 100-230/50-60                                                                              |
| Power output                            | W     | 250                                                                                        |
| Fuse (time-delayed)                     | A     | 3.15 / 6.3                                                                                 |
| Axes                                    | -     | 5                                                                                          |
| Accuracy                                | µm    | < 10                                                                                       |
| Allowable temperature range (operation) | °C    | 18 - 30                                                                                    |
| Spindle torque                          | Ncm   | 9.2                                                                                        |
| Spindle speed                           | rpm   | 100000                                                                                     |
| Collet diameter                         | mm    | 3                                                                                          |
| Sound pressure level                    | dB(A) | 60                                                                                         |
| Interfaces                              | -     | USB / Ethernet / slot for SD card / controlling of supply unit / controlling of extraction |
| Max. pump flow rate                     | l/min | 3.8                                                                                        |

Tab. 1

|                                     | Unit              | Value                                                                                                      |
|-------------------------------------|-------------------|------------------------------------------------------------------------------------------------------------|
| Compressed-air connection           | -                 | Dry, clean compressed air                                                                                  |
| _ Pressure                          | bar               | 6                                                                                                          |
| _ Flow rate                         | l/min             | 50                                                                                                         |
| Air purity                          |                   |                                                                                                            |
| _ Solid contaminants                | -                 | Class 3; better than 5 µm for solids                                                                       |
| _ Water content                     | -                 | Class 4; Max. pressure dew point +3 °C                                                                     |
| _ Total oil content                 | -                 | Class 2; max. oil content 0.1 mg/m <sup>3</sup>                                                            |
| Extraction                          | -                 | Suitable for zirconium oxide dust particles, HEPA microfilter (97.97 %), filter class H12, dust category M |
| _ Diameter of extraction connection | mm                | 38                                                                                                         |
| _ Suction capacity                  | m <sup>3</sup> /h | 240                                                                                                        |

Tab. 1

## Technical data of the supply unit

|                              | Unit | Value           |
|------------------------------|------|-----------------|
| Art. No.                     | –    | 178641          |
| Dimensions<br>(D × W × H)    | mm   | 500 × 200 × 400 |
| Weight                       | kg   | 5.2             |
| Capacity of C/L<br>container | l    | 20              |
| Mesh size of filter          | µm   | 50              |
| Interfaces                   | –    | Milling machine |

Tab. 2

## Accessories/Spare parts

| Art. No.   | Designation                                 |
|------------|---------------------------------------------|
| 178640     | Set-up cart with supply unit (package)      |
| 178641     | Supply unit                                 |
| 178642     | Set-up cart                                 |
| 179210     | Service unit                                |
| 179218     | Service set for spindle                     |
| 179230     | Collet, 3 mm                                |
| 181360     | Blank holder 98mm                           |
| 179253     | Swarf sieve                                 |
| 179254     | Extraction hood                             |
| 179255     | Splash guard                                |
| 179256     | Filter bag                                  |
| 179257     | Extraction hose with brush                  |
| 178650-NTR | Coolant/lubricant concentrate               |
| 178652     | pH-value test strips, 100 pce.              |
| 178653     | Nitrite test, 100 pce. test strips          |
| 1156901    | Fuse, milling machine<br>(T 3.15 A / 250 V) |

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**Amann Girrbach AG**
**Amann Girrbach GmbH**