
# TENAX Fiber Trans<sup>®</sup>

Esthetic Post System

Instructions for Use  
One-Office-Visit Technique

CE 0344  
 COLTENE

## TENAX Fiber Trans®

<b>EN</b>	Instructions For Use	4
<b>ES</b>	Instrucciones de uso	8
<b>PT</b>	Instruções de utilização	12
<b>ZH</b>	使用说明	16
<b>RU</b>	Инструкция по применению	19
<b>PL</b>	Instrukcja użycia	23
<b>CS</b>	Návod k použití	27
<b>SK</b>	Návod na použitie	31
<b>HU</b>	Használati útmutató	35
<b>SL</b>	Navodila za uporabo	39
<b>HR</b>	Upute za upotrebu	43
<b>RO</b>	Instrucțiuni de utilizare	47
<b>BG</b>	Инструкции за употреба	51

**TENAX** Fiber Trans®

### One-Office-Visit Technique

**Intended Use:** Fiber posts are intended to be cemented into the root canal of a tooth to stabilize and support a restoration.

 <b>CAUTION!</b>	Federal (United States) law restricts this device to sale by or on the order of a healthcare practitioner.
 <b>WARNING!</b>	Failure to properly follow the directions in this IFU may result in patient injury or compromise the integrity of the restoration.
 <b>NON STERILE</b>	Posts are delivered non-sterile and are to be sterilized prior to use.
 <b>WARNING!</b>	Material is composed of a methacrylate-based glass fiber reinforced resin. Device should not be used in patients with known material sensitivity.

Note: Posts are single use devices and should never be salvaged for re-use. Structural integrity of the post could be compromised including microscopic damage which could lead to post fracture.

For illustrations see page 55

**We strongly recommend the use of Hygenic® Rubber Dam with this procedure.**

1. After endodontic therapy has been completed, prepare the tooth for restoration as if pulp were intact. The preparation should include at least 1.5 mm of sound tooth structure all around the circumference of the preparation, apical to the core, for desired ferrule effect (**Fig. 01**). \*If that amount of tooth structure is not available, it should be created by periodontal crown lengthening and/or orthodontic extrusion.\*
2. Use a radiograph to determine the appropriate diameter and depth of the planned post space preparation. Keep in mind that sufficient root wall thickness must be maintained (1 mm minimum)\* to prevent perforation or weakening of root walls, and at least 4 - 5 mm of gutta-percha must be left intact to protect the apical seal.
3. Remove gutta-percha to preplanned depth with a Gates-Glidden drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended. If the canal has been obturated with a silver cone, remove it and reseal (obturate) with Hygenic gutta-percha before preparing the post space.
4. To begin creating the tapered post space, select the TENAX Fiber Trans Drill diameter which corresponds to the last Gates-Glidden drill used in Step 3 (See Approximate Diameter Comparison chart **Fig. A**). If using the TENAX Fiber Trans post as a parallel post, see NOTE: Using TENAX Fiber Trans as a Parallel Post below.

TENAX Fiber Trans Drills can be used either manually with the Universal Hand Driver or with a slow-speed contra-angle (750 RPM - 1,000 RPM). When a contra-angle is used, the drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the risk of the drill seizing in the post space. Remove all weakened or unsupported tooth structure. Remove any tooth debris from canal by irrigating post space with water.

5. Sequentially step up to the next larger TENAX Fiber Trans Drill until the preplanned diameter and depth are achieved.
6. Use a cylindrical diamond or carbide bur to prepare an anti-rotational box (**Fig. 02**).
7. Select the TENAX Fiber Trans Post that corresponds to the last drill used in the preparation of the post space. Insert post into post space (**Fig. 03**). Carefully check for occlusal clearance. Remove post; shorten as necessary from the non-tapered (occlusal) end (as clinical judgment dictates) with a Carborundum separating disk. Wear full face safety mask when shortening post. (Wet the post with water before cutting and rotate the post when cutting to ensure an even cut through each fiber.)
8. After try-in and before cementing, gently clean post with an alcohol wipe.
9. (Optional) Brush silane onto the post and allow to dry.
10. Condition the canal and/or apply primer/adhesive as per cement

manufacturer's instructions (we recommend ParaCore™).

11. Use the cement components according to manufacturer's instructions and apply to post surface. Use a Lentulo Spiral to coat the walls of the post space with cement, if setting time of cement will allow (**Fig. 04**).
12. Seat post immediately, slowly inserting the post to full depth, allowing excess cement to vent (**Fig. 05**). Apply pressure for about 60 seconds.
13. In the event light-curing is required, wipe away any excess cement before activating the light and follow manufacturer's instructions. Place curing light probe directly over the post, applying gentle pressure with the curing light tip. Light-cure for a minimum of 60 seconds (depending on volume and depth). If you are using an LED light, make sure your light-cured cement contains the Camphorquinone (CQ) Photoinitiator. If you have any questions, contact the cement manufacturer.
14. Complete the core and fabricate the final restoration (**Fig. 06**).

**NOTE: Using TENAX Fiber Trans as a Parallel**

**Post:** Follow steps 1-14, substituting ParaPost Drills (non-tapered, parallel drills) for the TENAX Fiber Trans Drills according to the Drill Application Comparison Chart (**Fig. B**). The TENAX Fiber Trans post will be placed with the non-tapered end in the apical portion of the canal, and shortened if necessary from the tapered end.

## Cleaning and Sterilization

**Posts:** Initiate cleaning by soaking the post(s) for 5 minutes in a solution consisting of an enzymatic cleaner, such as Terg-a-Zyme and de-ionized water at approximately 20°C (68°F). Follow the agent manufacturer's instructions.

Clean the post(s) by brushing away from the post while keeping it immersed in an enzymatic cleaning/disinfecting agent validated for use on dental instruments. Special care should be taken to clean crevices and other hard to reach areas. Thoroughly rinse with tap water. Dry using a paper towel.

Additional cleaning can be performed in an ultrasonic cleaner for 5 minutes. Remove and thoroughly rinse with tap water; dry with a paper towel.

Place post(s) into a small, clean covered glass container.

**Autoclave Sterilization:** ParaPost Fiber Trans may be dynamic-air-removal steam sterilized (4 minutes at 132°C followed by drying for 20 minutes). Follow manufacturer's recommendations. Posts must be sealed in a sterilization pouch before use in autoclave. Alternatively, sterilize in a pre-vacuum chamber at 132 °C for 4 minutes, 134 °C for 3 minutes or 134 °C for 6 minutes.

**Drills:** Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush; thoroughly rinse and dry. Place drills, with a sterilization indicator, into an autoclave bag and seal. Steam sterilize at 132°C for 12 minutes. Sterilizing bags should not touch the walls

of the autoclave. After use, to decontaminate, soak in a 70% isopropyl alcohol solution for 10 minutes. Once the sterilization monitor indicates that sterilization has been achieved, place autoclave bag, with the drills inside, in storage until ready for use.

**Universal Hand Driver:** Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush. Sterilize by following instructions for drills above.

**Warning Notice: DO NOT BEND POSTS!**

## Post Removal Technique

1. Take a radiograph of the tooth that requires post removal.
2. Carefully estimate length and diameter of TENAX Fiber Trans post to be removed. If possible consult patient's record.
3. Remove any restorative or core material which may be in the way to gain direct straight-line access to the post.
4. Wear full face safety mask when removing post from tooth. Indent the center of the post using a #1/2 round bur at high speed.
5. Carefully drill an initial channel into the center of the TENAX Fiber Trans post using a Kodex twist drill (cat. # K95 or K97) at slow speed.
6. Carefully drill to full depth through the center of the post using an end cutting Tenax Starter Drill (cat # TEDC1,

Ø .036"/0.90 mm) or Tenax Drill (cat # TEA11, Ø. 043"/1.1 mm). Check radiographically at each stage.

7. Gradually work up to original diameter using standard TENAX Fiber Trans drills.
8. Re-treat endodontically, if necessary.
9. Re-restore with new post/core and restoration.

Always verify steps via radiograph.

\*References and/or literature available upon request.

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




Fig. A Approximate Diameter Comparisons

Tenax Fiber Trans Drill		ISO D3 Instrument Measurement	Peeso Measurement mm	Gates-Glidden Measurement mm
Size (mm)	Grooves			
1.10	1	#60 0.92	#2 0.85	#3 0.90
1.30	3	#90 1.22	#4 1.22	#4 1.10 #5 1.30
1.50	5	#110 1.42	#5 1.40 #6 1.50	#6 1.50

**Fig. B Drill Application Comparisons**

<b>Tapered Application</b>		<b>Parallel Application</b>	
Tenax Fiber Trans Drill		ParaPost Drill	
Size	Drill Identification	Size	Drill Identification
1.1 mm	1 Groove	#4.5	Blue
1.3 mm	3 Grooves	#5	Red
1.5 mm	5 Grooves	#6	Black

## TENAX Fiber Trans Order Information

### TFT1 Tenax Fiber Trans One-Office-Visit Introductory Kit

Contents	5 Tenax Fiber Trans Posts	1.1 mm
	1 Tenax Fiber Trans Drill	1.1 mm
	5 Tenax Fiber Trans Posts	1.3 mm
	1 Tenax Fiber Trans Drill	1.3 mm
	5 Tenax Fiber Trans Posts	1.5 mm
	1 Tenax Fiber Trans Posts	1.5 mm
	1 Universal Hand Driver	

### REFILLS

5 Tenax Fiber Trans Posts per refill

TFT11	1.1 mm
TFT13	1.3 mm
TFT15	1.5mm

### POST REMOVAL:

TEDC1 Starter Drill and Universal Hand Driver

K95 Kodex Twist Drill

# TENAX Fiber Trans®

## Esthetic Post System

Coltène/Whaledent Inc.

235 Ascot Parkway

Cuyahoga Falls, OH 44223 / USA

Tel. USA & Canada +1 800 221 3046

+1 330 916 8800

Fax

+1 330 916 8800

info.us.coltene.com

Coltène Whaledent AG

Feldwiesenstrasse 20

9450 Altstätten / Switzerland

Tel. +41(0)91 757 53 00

Fax +41(0)71 757 53 01

info.ch.coltene.com

Coltène/Whaledent GmbH + Co. KG 

Raiffeisenstrasse 30

89129 Langenau / Germany

Tel. +49 (0) 7345 805 0

Fax +49 (0) 7345 805 201

info.de@coltene.com

www.coltene.com

Made in USA

 0344  

P/N: 30004079D 07.20.2018