ProTaper Next®

EN	Instructions For Use	2
FR	Instructions d'utilisation	7
IT	Istruzioni per l'uso	13
DE	Gebrauchsanweisung	19
ES	Instrucciones de uso	25
PT	Instruções de Utilização	30
NL	Gebruiksaanwijzing	35
SV	Användarinstruktioner	40
DA	Brugsanvisning	45
NO	Brukerinstruksjoner	50
LT	Naudojimo instrukcija	55
LV	Lietošanas norādījumi	60
ET	Kasutusjuhised	65
CS	Návod k použití	70
SK	Návod na použitie	75
HU	Használati utasítás	80
RO	Instrucțiuni de utilizare	85
SR	Uputstva za upotrebu	90
BS	Uputstva za upotrebu	95
CNR	Uputstva za upotrebu	100
SQ	Udhëzimet për përdorim	105
EL	Οδηγίες χρήσης	110
TR	Kullanım Kılavuzu	116
PT-BR	Instruções de uso	121
KO	사용 지침	126

ProTaper Next®

EN

*FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
SINGLE USE - REPROCESSING PROHIBITED
STERILE – STERILIZED BY RADIATION*

INSTRUCTIONS FOR USE

SINGLE USE STERILE ENDODONTIC INSTRUMENTS

0) INDICATIONS FOR USE

- The products are indicated for treatment of endodontic disease.
- Intended purpose: Engine-driven instrument intended for root canal preparation (shaping and debridement of the root canal).
- Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as Endo specialists (Endodontist) and Dental Assistants.
- Instruments shall be used in combination with an endodontic motor (manually in severe curvatures).

1) CONTRAINDICATIONS

None.

2) CONTENT, COMPOSITION AND COMPATIBLE DEVICES

ProTaper Next® instruments are made of three main components: the working part made of a nickel-titanium alloy (M-Wire®) (except for instrument XA which is made of standard NiTi), a silicone stopper and a plated brass shank with colored rings (except for instrument XA).

The ProTaper Next® range comprises following endodontic instruments:

Code	Product	ISO XXX YYYv XXX-ISO size YYY-taper size v-variable taper	Available lengths	Packaging and content
A 0803	ProTaper Next®	X1: ISO 017 004v X2: ISO 025 006v X3: ISO 030 007v X4: ISO 040 006v X5: ISO 050 006v	21/25/31 mm	Blister with 3 instruments (sterile)
	ProTaper Next® Accessory Opener File	XA: ISO 019 035	19 mm	Blister with 3 instruments (sterile)
	ProTaper Next® Assortment	X1 X2 X3	21/25/31 mm	Blister with 1 instrument of each ISO (sterile)
A 0800	ProTaper Next®	XA: ISO 019 035 X1: ISO 017 004v X2: ISO 025 006v X3: ISO 030 007v X4: ISO 040 006v X5: ISO 050 006v	21/25/31 mm	Blister with 6 instruments (sterile)
	ProTaper Next® Accessory Opener File	XA: ISO 019 035	19 mm	Blister with 3 instruments (sterile)
	ProTaper Next® Assortment	X1 X2 X3	21/25/31 mm	Blister with 2 instruments of each ISO (sterile)
A 0802	ProTaper Next® DEMO	X1 X2 X3	25 mm	Blister with 1 instrument of each ISO (sterile)

- Torque control devices ensure optimal usage.
- Use with endodontic motors in a constant rotation at a speed of 300 rpm.
- Set torque at:

File Size	Torque [Ncm]	Speed [rpm]
ProTaper Next® X1 ProTaper Next® X2 ProTaper Next® X3 ProTaper Next® X4 ProTaper Next® X5	2-3	300
ProTaper Next® XA	3-4	

- Lubricants such as NaOCl, EDTA, ProLube, Glyde™ shall be used.
- The use of radiographs in combination with an apex locator and a tool for adjusting the silicon stopper to the correct working length is the appropriate method of working length determination.

3) WARNINGS

- Strictly follow this Instruction for Use to minimize following risks to the device, the patient and/or the user:

- Breakage of instrument.
- Cross-contamination.
- Heat generation due to insufficient lubrication and irrigation.
- Swallowing of working part of the instrument.
- Toxic or allergic reactions caused by processing residues.
- Similar to all machine-driven root canal instruments, ProTaper Next® instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation. In this case, pre-curved hand files should be used in the apical region.

4) PRECAUTIONS

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Before using any instrument, make sure it is well connected to the contra-angle head.
- Check instrument and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Irrigate abundantly and frequently the canal throughout the procedure and after each instrument used (according to good dental practices).
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).
 - Use manual instead of rotary files.
 - Use small size, flexible or/and NiTi files (this will help avoid canal transportation).

- Visually inspect the working part for all the defects listed in the former paragraph during use (i.e after each wave).
- Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

5) ADVERSE REACTIONS

In the present technical state, no adverse reaction has been reported so far.

6) STORAGE CONDITIONS

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

7) STEP BY STEP INSTRUCTIONS

- 1) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
- 2) Create straight-line access to the canal orifice(s) with emphasis on flaring, flattening, and finishing the internal axial walls.
- 3) Irrigate abundantly and frequently after each instrument used (according to good dental practices).
- 4) Scout the whole root canal with small-sized hand files, determine working length by using an apex locator in combination with radiographs, irrigate and confirm patency.
- 5) If preflaring of the orifice is necessary, use the ProTaper Next® XA (Accessory) file to improve radicular access by removing triangles of dentin, relocating the coronal aspect of a canal away from an external root concavity and creating more shape as desired.
- 6) Always irrigate and if necessary, expand the glide path using small-sized hand files or dedicated mechanical glide path files.
- 7) In the presence of NaOCl, brush and follow along the glide path with the ProTaper Next® X1 (ISO 017 004), in one or more passes, alternatively with small-sized hand files if necessary, until full working length is passively reached.
- 8) During protocol of use, irrigate, recapitulate with a small-sized hand file after each sequential ProTaper Next® instrument, then re-irrigate.
- 9) Use ProTaper Next® X2 (ISO 025 006), exactly as described for ProTaper Next® X1, until full working length is passively reached.
- 10) Inspect the apical flutes of the ProTaper Next® X2: If they are loaded with dentin, then the shape is finished, the correspondingly sized gutta-percha master cone or size verifier may be fitted, and the canal is ready for disinfection.
- 11) Alternatively, gauge the foramen with an ISO 025 hand file:
 - If the instrument is snug at length, preparation is finished and the canal is ready to be obturated.
 - If the instrument is loose at length, continue shaping with ProTaper Next® X3 (ISO 030 007) and, when necessary, proceed with ProTaper Next® X4 (ISO 040 006) or ProTaper Next® X5 (ISO 050 006), gauging after each instrument respectively with a hand file of the same size. Use the same shaping technique up to working length.

12) Use dedicated paper points to dry and dedicated gutta-percha points to obturate the root canals.

8) HYGIENE, DISINFECTION, CLEANING AND STERILIZATION

- ProTaper Next® instruments are single use devices and cannot be reused. Re-use can increase the risk of cross contamination or breakage.
- Disinfection, cleaning and sterilization is not applicable.
- Products shall be disposed according to local regulations for the safe disposal of sharp and contaminated devices.

9) ADDITIONAL INFORMATION

- Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
- Sterility cannot be guaranteed if packaging is open, damaged or wet.
- To get a free printed copy of IFU please see section “Get a print copy of IFU” on website <https://www.dentsplysirona.com/en/service-contact/ifu.html>.
- Explanation of non-harmonized symbols for IFUs and labels, see IFU Symbols (<https://www.dentsplysirona.com/en/service-contact/ifu/ifu-symbols.html>).

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DENTSPLY DeTrey GmbH