GLUMA 2 Bond

(GB)	Instructions for use	3	(TR)	Kullanma talimatı	77
(ES)	Instrucciones de uso	10	(RU)	Инструкция по применению ..	84
(GR)	Οδηγίες χρήσης	18	(AE)	تعليمات الاستخدام	92
(CZ)	Návod k použití	26	(KR)	사용 설명서	98
(HU)	Használati utasítás	33	(SK)	Návod na použitie	105
(LV)	Lietošanas instrukcija	41	(BG)	Инструкции за употреба	113
(LT)	Naudojimo instrukcija	48	(UA)	Інструкція із застосування	121
(PL)	Instrukcja obsługi	55	(EE)	Kasutusjuhend	129
(HR)	Upute za uporabu	63	(SI)	Navodila za uporabo	136
(RO)	Instrucțiuni de utilizare	70			

GLUMA® 2 Bond

Instructions for use

Medical Device – for use by dental health care professionals only.

Intended purpose

GLUMA 2 Bond is a light-cured single-component adhesive for use in adhesive restorative dentistry. GLUMA 2 Bond has been developed for the adhesive bonding of plastic filling materials (e.g. composites, compomer) to dental hard substances and laboratoryfabricated restorations (e.g. ceramics). With GLUMA 2 Bond, priming, bonding and desensitizing can be done in one step. GLUMA 2 Bond is a light-activated adhesive in an ethanolic solution. Before using GLUMA 2 Bond, the dental hard substance is conditioned with an etching gel (etch&rinse).

Indications

- Adhesive fixing of direct composite and compomer restorations
- Adhesive fixing of indirect laboratory-fabricated ceramic and composite restorations (inlays, onlays, veneers and crowns)
- Treatment of hypersensitive tooth regions

Contraindications

In case of known or suspected allergy to components of the product, its use is contraindicated. Do not use this product in case of known or suspected allergies against (meth)acrylate compounds or glutaraldehyde. Do not bring GLUMA 2 Bond into direct contact with opened pulp tissue.

Patient target group

This dental material is suitable for patients requiring dental treatment for the following indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.

Intended user

For use by dental health care professionals only.

Clinical benefit

Kulzer medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthetically.

Composition

GLUMA 2 Bond contains: UDMA, ethanol, HEMA, water, methacrylated polycarboxylic acid, 4-META, glutaraldehyde

Potential side effects

Hypersensitivities to the product or its components cannot be excluded in individual cases.

Product features

Available in a 4 ml bottle

Application**1. Adhesive fixing of direct composite and compomer restorations****1.1 Preparation**

- Clean tooth with a toothpaste free from oil and fluoride.
- Prepare cavity as specified by the rules of adhesive filling treatment.
- Rinse prepared site with water and blow-dry. Avoid surface contamination.
- For deep cavities, we recommend protecting the dentine near the pulp with an appropriate base filling (e.g. calcium hydroxide preparation and glass ionomer cement base filling).
- Cover the dentine near the pulp only. Leave the remainder of the cavity free for adhesion of GLUMA 2 Bond.

1.2 Conditioning of the tooth surface

Before applying GLUMA 2Bond, the prepared cavity must be pretreated with a suitable etching gel according to the manufacturer's instructions (e.g. GLUMA Etch 35 Gel from Kulzer) and thoroughly rinsed with water and dried. Overdrying of the dentin should be avoided; the dentin surface should be slightly damp.

1.3 Application of GLUMA 2Bond

- Dispense GLUMA 2Bond in a mixing well, extract it with an applicator tip or a soft disposable brush and apply to the entire surface of the cavity.
- Allow GLUMA 2Bond to soak in for 15 seconds.
- GLUMA 2Bond contains ethanol as a solvent. Carefully evaporate solvent and residual moisture under a gentle air stream until no movement of liquid can be detected.

An excessively strong air flow at the start of blow-drying will dilute the bonding and may result in inadequate adhesion.

- The surface must be visibly shiny after application of GLUMA 2Bond and after evaporation of the solvent. Make certain that the complete surface of the cavity is fully covered. If the cavity does not appear to be completely shiny all over, apply GLUMA 2Bond a second time as described above. Make sure that no excess remains in corners of the cavity. Cure GLUMA 2Bond for 20 seconds using a suitable halogen or LED polymerisation lamp (wavelength 440–470 nm, light intensity at least 600 mW/cm²). If the light power is too low there will be insufficient adhesion.

Insufficient lamp intensity will result in inadequate adhesion. Polymerisation lamps should be tested with reliable testing devices at regular intervals. The light emitting window should be placed as close as possible to the plastic surface during polymerisation.

- Then immediately apply and process the filling material in the cavity as directed by the manufacturer.

2. Adhesive fixing of indirect laboratory-fabricated ceramic and composite restorations (inlays, onlays, veneers and crowns)

2.1 Preparation of the restoration

- Prepare restoration surface as directed by the manufacturer's instructions for use.

2.2 Conditioning of the tooth surface

Proceed as described in 1.2

2.3 Application of GLUMA 2Bond

Make sure that no excess material is left in the corners of the cavity and that the GLUMA 2Bond coating is carefully blow-dried before polymerisation.

Otherwise proceed as described in 1.

- Polymerise GLUMA 2Bond for 20 s with a halogen or LED polymerisation lamp (see instructions above).
- Place composite cement (self-curing or dual-curing) according to the manufacturer's instructions for use.

3. Treatment of hypersensitive tooth regions and cleaning the tooth

- Carefully rinse the surface with water and blow-dry.

3.1 Conditioning hypersensitive tooth regions

Proceed as described in 1.

3.2 Application of GLUMA 2Bond

Proceed as described in 1.

- Polymerise GLUMA 2Bond for 20 s with a halogen or LED polymerisation lamp (see instructions above).
- Carefully remove the oxygen inhibition layer with an alcohol-soaked pellet.

If GLUMA 2Bond does not sufficiently desensitise, repeat the process as described above.

General information

Please pay attention to the following safety instructions and all processing instructions in other sections. This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself.

- Visually inspect for damage before use. Damaged products must not be used.
- Materials containing eugenol may impair the polymerisation of GLUMA 2Bond.
- Ensure product is at room temperature prior to application.

Warnings and safety instructions

Highly flammable liquid and vapour. May cause skin and eye irritation. Avoid contact with skin. Skin contact may cause sensitization. Wear suitable protective gloves. In case of contact with skin, immediately wash with plenty of water and soap. Avoid contact with eyes. In case of contact with eyes, rinse with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If symptoms persist, consult an eye specialist. Toxic to aquatic life with long lasting effects.

- GLUMA 2Bond contains ethanol. Ethanol is very flammable. Do not inhale vapour.

- GLUMA 2Bond contains glutaraldehyde. If it is inhaled or swallowed, it may damage health. Seek medical attention immediately if swallowed. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources - no smoking. May cause respiratory irritation. Avoid breathing dust/fume/gas/mist/vapours/spray. Avoid static charging. Keep away from all sources of ignition.
- Use only in sufficiently well-ventilated rooms. Do not inhale vapours.
- Do not apply adhesive directly to the exposed pulp.
- Eye protection is recommended for the patient. The use of a rubber dam is strongly recommended.
- Close the GLUMA 2Bond bottle carefully after each use.
- To prevent cross-contamination, do not reuse bottles after contamination with bodily fluids.
- Saliva or blood in the field of application during treatment may lead to failure. Use rubber dam (e.g. Ivory®) or other suitable means of isolation.

Special advice

For SSCP please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact sscp-service@kulzer-dental.com

To find the SSCP in the EUDAMED use the UDI code attached at the packaging of the product. According to EU Medical Device Regulation users/ patients are obliged to report serious incidents with a medical device to the manufacturer and to the competent authority of the country, the incident occurred.

Storage conditions

Do not use the material after the expiration date. After opening, use by expiry date. Store at temperatures of -10 to 25°C (14 to 77°F). The product can be stored in refrigerator. Use the adhesive at room temperature only. Do not expose to direct sunlight. Store in a sufficiently well ventilated dry place. The shelf life will be reduced by unsuitable storage conditions and the function of the product may be affected. Close the container tightly immediately after use.

Disposal information

Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system. European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

Complaint notice

Please quote batch number and article number or UDI in all correspondence about the product. Article number **REF** and batch number **LOT** or UDI: See note on bottle/packaging.

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Kulzer GmbH