Prior to use, carefully read the instructions for use. **EN**

Saliva-Check BUFFER

In Vitro Test for checking the Quality, pH and Buffering Capacity of Saliva

For use only by a dental professional in the recommended indications.

INDICATIONS FOR USE

To check the properties of resting and stimulated saliva.

CONTRAINDICATIONS

In vivo diagnostic use.

DIRECTIONS FOR USE

Prior to any visit where a saliva diagnostic procedure is planned, instruct the patient not to smoke, consume food or drink, brush the teeth or use a mouth wash for at least one hour prior to the scheduled appointment time.




I. TESTING OF RESTING SALIVA

TEST 1 – Visual inspection of level of hydration Visually assess the lower lip labial gland secretion. Evert the lower lip, gently blot the labial mucosa with a small piece of gauze and observe the mucosa under good light. Droplets of saliva will form at the orifices of the minor glands. Assess the time for visible production of saliva as follows¹:


Greater than 60 seconds : resting flow Low 

Less than 60 seconds : resting flow Normal 

TEST 2 – Saliva consistency Visually assess the resting salivary consistency in the oral cavity¹.

Sticky frothy saliva residues: Increased viscosity 
Frothy bubbly saliva: Increased viscosity 
Watery clear saliva: Normal viscosity 

TEST 3 – pH measurement Instruct the patient to expectorate any pooled saliva into the collection cup. Take a pH test strip, place this into the sample of resting saliva for 10 seconds, and then check the colour of the strip. This should be compared with the testing chart available in the package.

II. TESTING OF STIMULATED SALIVA

TEST 4 – Saliva quantity Instruct the patient to chew the piece of wax to stimulate salivary flow. After 30 seconds, let the patient expectorate into the spittoon. Continue chewing for a further 5 minutes, collecting all the saliva into the collection cup at regular intervals.

The quantity of saliva can be measured by checking the mL markings on the side of the cup. Quantity of saliva at 5 minutes²:

< 3.5 mL Very low 
Between 5.0 – 3.5 mL Low 
> 5.0 mL Normal 

Note: Normal stimulated saliva flow rate may vary between 1 mL/min – 1.6 mL/min.


TEST 5 – Buffering capacity

- Remove a Buffer test strip from the foil package and place onto an absorbent tissue with the test side up.
- Using a pipette, draw sufficient saliva from the collection cup and dispense one drop onto each of the 3 test pads. Immediately turn the strip 90° to soak up excess saliva on the absorbent tissue. This will prevent the excess saliva from swelling on the test pad and possibly affecting the accuracy of the test result.
- The test pads will begin to change colour immediately and after 2 minutes the final result can be calculated by adding the points according to the final colour of each pad. See conversion table and examples underneath.

Conversion table

Test pad colour at 2 minutes	
Green	4 points
Green / Blue	3 points*
Blue	2 points
Red / Blue	1 point *
Red	0 points


Examples

*Where a colour combination provides an unclear result, use intermediate scores.

Interpreting the result.

Combined total Buffering ability of Saliva

0-5 Very low 
6-9 Low 
10-12 Normal / High 

STORAGE

Recommended for optimal performance, store components at room temperatures (15-30°C).

PACKAGE

Contents	Quantity
In vitro pH test strips	20
Saliva collection cup	20
Wax for saliva stimulation	20
Saliva dispensing pipette	20
Buffer test strip	20

CAUTION

- In warm ambient temperatures, store the wax in a refrigerator.
- After opening the wax and Buffer test strip, use immediately.
- All components used for each test are for single use only. Do not reuse.
- Dispose of according to local regulations.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> or <http://The Americas.gceurope.com>. They can also be obtained from your supplier.

*1 Information by kind permission of Prof. L Walsh©2002
*2 Lund University Faculty of Odontology Department of Cariology, Sweden

Undesired effects- Reporting: If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

UK Responsible Person
GC UNITED KINGDOM Ltd.
Coopers Court Newport Pagnell, Buckinghamshire, MK16 8JS
United Kingdom

Last revised: 04/2022

Rx Only

DISTRIBUTED BY
GC CORPORATION
7-1-1 Hasegawa-cho, Hasegawa-ku, Tokyo 174-8585, Japan
EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00
GC AMERICA INC.
3737 West 127th Street, Aliso Viejo, CA 92683 U.S.A.
TEL: +1-714-557-0500
www.gcamerica.com
GC SOUTH AMERICA
Rua Heliodora, 399, Santana - São Paulo, SP, BRASIL
CNPJ: 02022-019 - TEL: +55-11-2925-9995
CNPJ: 08.279.999/0001-61
RESP. TEC: Etick de Lima - CRCSP 100.886

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01, Tampines Plaza 2, Singapore 529541
TEL: +65 6546 7168
GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Bankmeadow, NSW 2019, Australia
TEL: +61 2 9301 8200
MANUFACTURED BY
NIMI CHEMICAL INDUSTRIAL CO., LTD.
4-13 Yokoyama-cho, Kiryu, Gunma, Japan
RESPONSIBLE MANUFACTURER IN CANADA:
GC AMERICA INC.
3737 West 127th Street, Aliso, CA 92683 U.S.A.
PRINTED IN JAPAN

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