

Cementation Solutions

# Ketac™ Cem Plus Automix

Resin Modified Glass Ionomer Cement



Enjoy easy clean-up  
within seconds

**3M** ESPE

# Ketac™ Cem Plus Automix – comfort plus performance



Now even easier to use with  
**NEW** tack light-cure clean-up feature.

Ketac™ Cem Plus Resin Modified Glass Ionomer (RMGI) Cement from 3M™ ESPE™ offers the perfect combination of comfortable handling and proven clinical performance. It's well suited for the cementation of PFM restorations and of restorations to implant abutments.

The convenient Automix delivery system makes it easier to use.

There are good reasons for preferring the automix syringe: thanks to this leading-edge device, Ketac Cem Plus offers exceptionally easy handling. And with its unique new material properties and the easy cleanup of excess, you will find that resin modified glass ionomer cementation is now more convenient than ever.



[www.dentaladvisor.com](http://www.dentaladvisor.com)

## Benefits

### **NEW** TACK-CURE FEATURE!

- **Faster, easier clean up**
- **No clogging**
- **Fluoride release**
- **Customer preferred automix delivery system**



## Indications

- **PFM crowns and bridges**
- **Metal crowns, bridges, inlays and onlays**
- Endodontic posts
- Orthodontic appliances
- All-zirconia or all-alumina strengthened core-ceramic restorations (such as Lava™, Procera®, AllCeram)
- **NEW!** Cementation to implant abutments



Ketac™ Cem Plus Cement is recommended for the cementation of restorations to implant abutments.

# Enjoy an easy clean-up with tack light cure option

Enjoy easy clean-up of excess cement within seconds.


You can speed up the gel state by light-curing excess Ketac™ Cem Plus Cement for about 5 seconds per surface with the Elipar™ S10 Curing Light from 3M ESPE or any other conventional light-curing device. Alternatively, you can wait for the excess to self-cure – about two minutes after seating in the mouth. Once excess Ketac Cem Plus cement reaches a gel state, it can be very easily dislodged in large sections from the margins.



Tack curing of Ketac™ Cem Plus Cement to remove excess easily.



Easy removal of excess cement in large pieces.




Compared to competitive products, Ketac™ Cem Plus Cement shows high shear bond strength values to enamel/dentine.

## Bond strength.

Ketac Cem Plus cement provides molecular bonding to the tooth surface. Compared to other leading cements, it delivers exceptional shear bond strength values – without a separate tooth conditioner. And there is virtually no post-operative sensitivity.

Source: 3M ESPE internal data 2012 (Notched-edge adhesion).



Ketac™ Cem Plus Cement offers sustained fluoride release.

## Fluoride release.

Ketac Cem Plus cement offers sustained fluoride release.

Source: Partially cited from: Ton T, Wang Y, Falsafi A, Oxman J. Fluoride Release of a Novel Automix Resin-Modified Glass-Ionomer Luting Cement. J Dent Res 91 (Spec Iss A) 1319, 2012. (www.dentalresearch.org).

## Ordering Information



Ketac™ Cem Plus Automix

Item No.	Product Information
	<b>Ketac™ Cem Plus Automix</b>
<b>3536</b>	<b>Value Pack</b> 3 Syringes à 8,5g, 36 Mixing Tips
<b>3536TK</b>	<b>Trial Kit</b> 1 Syringe à 8,5g, 12 Mixing Tips
<b>3541</b>	<b>Mixing Tips Refill</b> 12 pieces

Also available in the Clicker™ Dispenser.



Ketac™ Cem Plus Clicker™

	<b>Ketac™ Cem Plus Clicker™</b>
<b>56931</b>	<b>Trial Kit</b> 1 Clicker™ Dispenser 11g, 1 mixing pad
<b>56930</b>	<b>Refill</b> 2 Clicker™ Dispenser 11g each, 1 mixing pad

**3M ESPE**

3M Deutschland GmbH  
Location Seefeld  
3M ESPE · ESPE Platz  
82229 Seefeld · Germany  
info3mespe@mmm.com  
www.3MESPE.com



3M, ESPE, Clicker, Elipar, Ketac and Lava are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada.

All other trademarks are owned by other companies.

© 3M, 2012. All rights reserved.

XX (0X.2012)