## Celtra® Ceram

Instructions for Use  
 Gebrauchsanweisung  
 Mode d'emploi  
 Istruzioni per l'uso  
 Instrucciones de uso  
 Instruções de Utilização  
 Οδηγίες χρήσης

Gebbruksaanwijzing  
 Bruksanvisning  
 Brugsanvisning  
 Bruksanvisning  
 Naudojimo instrukcija  
 Lietošanas instrukcijas  
 Kasutusjuhised

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## Celtra® Ceram

### Introduction

Thank you for choosing Celtra® Ceram veneering porcelain. For more than 100 years, Dentsply Sirona has been shaping the dental industry with innovative products designed to establish a new standard in esthetics and quality in fabricating all-ceramic dental prostheses. We are pleased to introduce this new product system with the highest level of performance our customers demand.

### Indications

Celtra Ceram is suitable for veneering all-ceramic frameworks.

### Compatibility

Celtra Ceram is designed exclusively for dental use only by trained professionals. Specifically, it is a low-fusing, leucite-reinforced feldspathic ceramic optimized for veneering and characterizing all ceramic frameworks (see below) in a dental laboratory.

With a coefficient of thermal expansion (CTE) of  $9.0 \times 10^{-6} \text{K}^{-1}$  (25-500 °C) and a firing temperature of 770 °C (1st dentin), Celtra Ceram is suitable for lithium disilicate frameworks. A firing temperature of 780 °C (1st dentin) is recommended for zirconia substrates.

- Celtra® Press Zirconia-reinforced Lithium Silicate (ZLS) frameworks:  
CTE  $9.7 \times 10^{-6} \text{K}^{-1}$  (25-500 °C)
- Lithium disilicate frameworks:  
CTE  $10.0 - 10.5 \times 10^{-6} \text{K}^{-1}$  (25-500 °C)
- Cercon ht zirconia framework:  
CTE  $10.5 \times 10^{-6} \text{K}^{-1}$  (25-500 °C)
- Cercon xt zirconia framework:  
CTE  $10.1 \times 10^{-6} \text{K}^{-1}$  (25-500 °C)
- Zirconia frameworks:  
CTE  $10.1 - 11.0 \times 10^{-6} \text{K}^{-1}$  (25-500 °C)

### Contraindications

Only the indications listed above are suitable. The following are not acceptable for use of Celtra Ceram:

- Celtra DUO blocks with a CTE of  $11.6 \times 10^{-6} \text{K}^{-1}$  (25-500 °C)
- Titanium or alloy frameworks with a similar CTE
- Alumina frameworks
- Any other veneering ceramic system
- Bruxism or other parafunctions
- Insufficient occlusal distance

### Warnings

Undesired side effects of these medical products are extremely rare with proper processing and use. In case of skin sensitization or rash, discontinue use and seek medical attention. Immunoreactions, such as allergies and/or local irritations (taste or oral mucosa), cannot completely be excluded on principle.

For patients with hypersensitivity to any of the ingredients, this medical device must not be used or only used under strict supervision of a treating physician/dentist. Similarly, known cross-reactions or interactions of this medical product with other materials existing in the mouth must be considered by the physician/dentist.

### Precautions

- Do not inhale abrasive dusts
- Do not ingest any of the materials (paste/powder/fluid)
- Frameworks that show evidence of cracking or surface voids must be discarded and not processed further; do not attempt to repair any cracks during the veneering process or additional firings
- Consider the following during firing to avoid potential for cracking:
  - Do not completely fill restorations with refractory putty
  - Use only ceramic pins, if possible
  - If wire pins are used, cover the wire with refractory putty but do not completely fill restoration
- Frameworks below minimum thickness are not recommended
- For minimum thickness recommendations for restorations, please consult table on page 7/8.

### Adverse reactions

No adverse reactions have been reported for Celtra Ceram. When working with these materials, make sure to comply with the Instructions for Use and the pertinent Safety Data Sheets (SDS).

### Technical data

- Coefficient of Thermal Expansion (CTE):  
 $9.0 \pm 0.5 \times 10^{-6} \text{K}^{-1}$  (25-500 °C),  $T_g = 520 \text{ °C}$
- Device classification: Dental ceramics, type 1, class 1 per ISO 6872
- Flexural strength: 108 MPa
- Chemical solubility: 28  $\mu\text{g}/\text{cm}^2$

## Handling, transport and storage

- Powders: Avoid direct exposure to sunlight and heat, protect from moisture, when not in use keep lid tightly closed.

## Symbols on product labels

-  Manufacturer
-  Batch Code
-  Reorder Number
-  Manufacture Date
-  Expiration Date
-  Consult Instructions for Use
-  Caution
-  Keep Dry
-  CE European Conformity Marking
-  Medical Device
-  U.S. Prescription Only

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

## Compatible materials

Multiple compatible materials and systems are recommended for use with the Celtra system and referenced herein. For more information related to these products, including ordering information, please visit [celtra-dentsplysirona.com](http://celtra-dentsplysirona.com).

## Compatible liquids

For best results, the following accessory liquids are recommended:

- Dentsply Sirona Modeling Liquid DU
- Dentsply Sirona Modeling Liquid U
- Ducera Liquid SD
- Dentsply Sirona Stain and Glaze Liquid

If longer working times are desired, the following liquid can be used as well:

- Dentsply Sirona Modeling Liquid E
- Ducera Liquid Form

## Isolation

- Ducera Isolating Fluid SEP
- Dentsply Sirona Die Release

## System Overview

Please remember when selecting shades that Celtra Press ingots and Cercon disks are made to the corresponding shade and will match the tooth shade exactly. The finished restoration will therefore correspond to the dentin shade and will require characterization only in the incisal area.

## Cut-Back Technique

### Enamels (E)

#### *E1-Extra Light, E2-Light, E3-Medium*

Enamels are optimized in translucency and color range to perfectly mimic natural effects in the incisal region

**6** **8** **9** **12** & **13**.

#### *Enamel Transparent (E4)*

Transparent powder without opalescence that can be used in a thin application to create depth and can be mixed with any other powders to increase translucency **14**.

#### *Enamel White (E5)*

Whitish effect powder to enhance occlusal cusps, palatal/lingual ridges or the incisal of anterior areas - can be diluted by using Enamel Opal Transparent EO4.

#### *Enamel Dark (E6)*

Enamel powder optimized in translucency and chroma designed exclusively for darker Shade Series dentin powders. See Shade Combination Table for guidance when Enamel Dark is recommended to maintain final shade.

### Enamels Opal (EO)

#### *EO1-Extra Light, EO2-Light, EO3-Medium*

Opalescent enamel inspired by nature to expand design options in the incisal area without affecting core color. Enamel Opal increases vitality without being too transparent (no greying effect) yet maintains a transparent visual character. A true multifunctional material to quickly achieve high esthetic results.

#### *Enamel Opal Transparent (EO4)*

A strong opalescent, almost transparent multi-functional powder with broad application within the system. EO4 can be used in pure form and/or mixed with all other powders to increase opalescence.

#### *Enamel Opal HT (EO5)*

A translucent enamel to enhance accents of yellow and orange opalescent effects. Especially formulated for highly translucent core materials, yet versatile to encourage creativity.

#### *Enamel Opal LT (EO6)*

An opalescent powder with higher value and higher opacity than EO5 intended for use with lighter and bleach shade substructures, yet similar to EO5 in versatility to promote creativity.

### Enamel Effect (EE)

Enamel Effects of varying hues that can be applied in incisal areas to enhance color depth and introduce natural features of adjacent dentition **2**, **3**, **4**, **7**, **10**, **11** & **13**.

#### *Enamel Effect Sunrise (EE1) / Enamel Effect Sunset (EE3)*

Opalescent effect powder can be used for yellow/orange (Sunrise) **7** & **13**, as well as orange/reddish (Sunset) accents in the enamel areas. They are well suited for increasing the chroma level at the 2nd or 3rd dentin firings. Sunset is mainly used for A-shades, whereas Sunrise is used mainly for B-shades. Both powders can be diluted using EO4 Enamel Opal Transparent to reduce the intensity of yellow and reddish colors.

#### *Enamel Effect Violet (EE2) / Enamel Effect Sky (EE5)*

Opalescent effect powders for discreet accent rod formations **2** & **11** and for creating strong illusions of depth and translucency in incisal areas. They can be diluted using EO4 Enamel Opal Transparent to reduce the intensity of the violet or sky.

#### *Enamel Effect Fog (EE4)*

Opalescent effect powder for grayish incisal areas **4** & **11** – can be diluted by using EO4 Enamel Opal Transparent to reduce the gray intensity.

#### *Enamel Effect Ivory (EE6)*

Whitish opalescent effect powder **3** & **10** for palatal/lingual ridges in the anterior area and to enhance occlusal cusps in the posterior region as well – can be diluted by using EO4 Enamel Opal Transparent to reduce the intensity of milky-ivory color.

## Layering Technique

### Dentins (D)

Dentins are optically balanced in hue, chroma, and value, with out-of-the-bottle shade accuracy and are available in a total of 46 shades: 4 BL1 – 4 bleach, 16 VITA® classical A1– D4 and 26 Shade Series – keyed to VITA® 3D-Master shade designations<sup>1</sup>.

### Power Dentins (PD)

Power dentins are highly chromatic and fluorescent ceramic masses for individual shade reproduction. They can be used in the cervical, palatal, or occlusal area to increase the chroma where required. This makes them particularly well-suited for veneering lithium disilicate frameworks. Depending on the desired effect and intensity, these masses can be used alone or in combination.

### Opaceous Dentins (OD)

Opaceous dentins resemble dentins in terms of hue and chroma, except that their opacity is increased by about 25%. Opaceous dentins can be used to control light reflection and translucency in the body of a ceramic restoration in situations where space is limited. They are primarily used for veneering zirconia frameworks. Depending on the desired effect and intensity, these masses can be used alone or in combination.

### Dentin Effects (DE)

Dentin Effects are intensive powders in a range of cool and warm hues which can be used as rod formations, accents and chromatic or value adjustments plus varied effects. EO4 Enamel Opal Transparent can be used to dilute and reduce the intensity.

### Dentin Gingiva (DG)

Gingiva powders in multiple shades mimicking soft tissue areas. Dentin Gingiva has same translucency as dentin and is fired during dentin applications at 770 °C, which is especially helpful in veneering implant-supported frameworks.

### Add-On Correction (C) & Add-On Gingiva (G)

Add-On Correction porcelain can be used for making final adjustments. Similarly, Add-On Gingiva can be used for final corrections in gingival areas. Both are low fusing ceramics with a firing temperature of 750 °C or 760 °C (see firing charts), and can be used either with or after glaze firing.

### Stains & Overglaze

Dentsply Sirona Universal Stains and Glaze has been especially developed with broad compatibility with a wide array of substructures and veneering systems, including Celtra Press, Cercon and Celtra Ceram restorations. More information about this universal system is available at [dentsplysirona.com/ceramics](https://dentsplysirona.com/ceramics).

**Note:** To more easily achieve full contour restorations in all 16 VITA® Classical shades, please see the Celtra and Cercon Staining Guides available at [dentsplysirona.com/ceramics](https://dentsplysirona.com/ceramics)

<sup>1</sup> VITA classical A1 – D4, and VITA® 3D-MASTER are registered trademarks of VITA Zahnfabrik H. Rauter GmbH & Co. KG.

## Shading recommendation – Celtra® Ceram shades keyed to Bleach and VITA Classical shades for Cercon® ht and Celtra® Press substrates/cores

Cercon ht core	Celtra Press ingot	Final Shade	Enamel (E)	Enamel Opal (EO)	Dentin	Opaceous Dentin	Power Dentin	Enamel Effect	Add-on Correction (C)
BL	BL1	BL1	E1-Extra light	EO6-LT	BL1	OD0	–	EE2-Violet, EE5-Sky	C1-Light
BL	BL1	BL2	E1-Extra light	EO6-LT	BL2	OD0 + OD1	–	EE2-Violet, EE5-Sky	C1-Light
BL	BL2	BL3	E1-Extra light	EO6-LT	BL3	OD0 + OD2	–	EE2-Violet, EE5-Sky	C1-Light
BL	BL2	BL4	E1-Extra light	EO6-LT	BL4	OD0 + OD2	–	EE2-Violet, EE5-Sky	C1-Light
A1	A1	A1	E1-Extra light	EO1-Extra light	A1	OD2	PD2	EE3-Sunset, EE2-Violet, EE5-Sky	C1-Light
A2	A2	A2	E1-Extra light	EO1-Extra light	A2	OD2 + OD5	PD3 + PD5	EE3-Sunset, EE2-Violet, EE5-Sky	C1-Light
A3	A3	A3	E2-Light	EO2-Light	A3	OD2 + OD5	PD3 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
A3.5	A3	A3.5	E3-Medium	EO2-Light	A3.5	OD5	PD4 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C3-Dark
A4	A3	A4	E3-Medium	EO3-Medium	A4	OD5 + OD6	PD5 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C3-Dark
B1	B1	B1	E1-Extra light	EO1-Extra light	B1	OD1	PD1	EE1-Sunrise, EE2-Violet, EE5-Sky	C1-Light
B2	B1	B2	E1-Extra light	EO1-Extra light	B2	OD2 + OD3	PD2 + PD3	EE1-Sunrise, EE2-Violet, EE5-Sky	C1-Light
B3	B3	B3	E1-Extra light	EO1-Extra light	B3	OD4 + OD5	PD2 + PD6	EE1-Sunrise, EE2-Violet, EE5-Sky	C2-Medium
B4	B3	B4	E2-Light	EO2-Light	B4	OD4 + OD5	PD4 + PD6	EE1-Sunrise, EE2-Violet, EE5-Sky	C3-Dark
C1	C1	C1	E1-Extra light	EO1-Extra light	C1	OD3	PD3	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
C2	C1	C2	E3-Medium	EO3-Medium	C2	OD3 + OD6	PD3 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
C3	C3	C3	E3-Medium	EO3-Medium	C3	OD3 + OD6	PD3 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C3-Dark
C4	A3	C4	E3-Medium	EO3-Medium	C4	OD6	PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C3-Dark
D2	D2	D2	E3-Medium	EO3-Medium	D2	OD3	PD2 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
D3	D3	D3	E3-Medium	EO3-Medium	D3	OD3 + OD5	PD3 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
D4	C1	D4	E3-Medium	EO3-Medium	D4	OD6	PD3 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium

**Opaceous dentin application for:** frameworks with low translucency, e. g. Cercon® base or Cercon® ht white

**Power dentin application for:** frameworks with high translucency, e. g. Cercon® ht, xt True Color Technology or Celtra® Press

## Shading recommendation – Celtra® Ceram shades keyed to VITA 3D-Master shades for Cercon® ht and Celtra® Press substrates/cores

Starting Core/ ingot shade		Final Shade	Enamel (E)	Enamel Opal (EO)	Dentin	Opaceous Dentin	Power Dentin	Enamel Effect	Add-on Correction (C)
Cetra Press (LT)	Cercon ht								
B1	B1	101 (1M1)	E1-Extra light	EO1-Extra light	101 (1M1)	OD1	PD1	EE3-Sunset, EE2-Violet, EE5-Sky	C1-Light
A1	A1	102 (1M2)	E1-Extra light	EO1-Extra light	102 (1M2)	OD2	PD1 + PD2	EE1-Sunrise, EE2-Violet, EE5-Sky	C1-Light
A1	A1	203 (2L1.5)	E1-Extra light	EO1-Extra light	203 (2L1.5)	OD2 + OD3	PD3	EE1-Sunrise, EE2-Violet, EE5-Sky	C1-Light
B3	B2	204 (2L2.5)	E1-Extra light	EO1-Extra light	204 (2L2.5)	OD2 + OD3	PD2 + PD3	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
A1	A1	205 (2M1)	E1-Extra light	EO1-Extra light	205 (2M1)	OD3	PD1 + PD3	EE3-Sunset, EE2-Violet, EE5-Sky	C1-Light
A2	A2	206 (2M2)	E1-Extra light	EO1-Extra light	206 (2M2)	OD2 + OD3	PD2 + PD3	EE1-Sunrise, EE2-Violet, EE5-Sky	C1-Light
A3	A3	207 (2M3)	E1-Extra light	EO1-Extra light	207 (2M3)	OD3 + OD4	PD2 + PD4	EE1-Sunrise, EE2-Violet, EE5-Sky	C2-Medium
A2	A2	208 (2R1.5)	E1-Extra light	EO1-Extra light	208 (2R1.5)	OD2 + OD3	PD2 + PD3	EE3-Sunset, EE2-Violet, EE5-Sky	C1-Light
A2	A2	209 (2R2.5)	E1-Extra light	EO1-Extra light	209 (2R2.5)	OD3 + OD5	PD1 + PD5	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
C3	C2	310 (3L1.5)	E2-Light	EO2-Light	310 (3L1.5)	OD1 + OD6	PD3	EE1-Sunrise, EE2-Violet, EE5-Sky	C2-Medium
B3	B2	311 (3L2.5)	E2-Light	EO2-Light	311 (3L2.5)	OD5	PD3 + PD4	EE1-Sunrise, EE2-Violet, EE5-Sky	C2-Medium
C1	C1	312 (3M1)	E1-Extra light	EO1-Extra light	312 (3M1)	OD1 + OD6	PD3	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
D3	D3	313 (3M2)	E2-Light	EO2-Light	313 (3M2)	OD3 + OD5	PD3 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
B3	B3	314 (3M3)	E3-Medium	EO3-Medium	314 (3M3)	OD4 + OD5	PD4 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
C1	C1	315 (3R1.5)	E2-Light	EO2-Light	315 (3R1.5)	OD3 + OD5	PD3 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C1-Light
D3	D3	316 (3R2.5)	E3-Medium	EO3-Medium	316 (3R2.5)	OD4 + OD5	PD3 + PD5	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
C3	C3	417 (4L1.5)	E6-Dark	2x EO3 + 1x EE4 Fog	417 (4L1.5)	OD3 + OD6	PD3 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C3-Dark
B3	B3	418 (4L2.5)	E3-Medium	EO3-Medium	418 (4L2.5)	OD5 + OD6	PD5 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C3-Dark
C3	C3	419 (4M1)	E6-Dark	2x EO3 + 1x EE4 Fog	419 (4M1)	OD3 + OD6	PD3 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
C3	C3	420 (4M2)	E6-Dark	2x EO3 + 1x EE4 Fog	420 (4M2)	OD5 + OD6	PD5 + OD6	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
A3	A3.5	421 (4M3)	E3-Medium	EO3-Medium	421 (4M3)	OD6 + PD5	PD5	EE3-Sunset, EE2-Violet, EE5-Sky	C2-Medium
C3	C3	422 (4R1.5)	E6-Dark	2x EO3 + 1x EE4 Fog	422 (4R1.5)	OD5 + OD6	PD5 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C3-Dark
C3	C3	423 (4R2.5)	E3-Medium	EO3-Medium	423 (4R2.5)	OD6 + DE9	PD5 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C3-Dark
C3	C3	524 (5M1)	E6-Dark	2x EO3 + 1x EE4 Fog	524 (5M1)	OD6	PD6 + OD6	EE3-Sunset, EE2-Violet, EE5-Sky	C3-Dark
C3	C4	525 (5M2)	E6-Dark	2x EO3 + 1x EE4 Fog	525 (5M2)	OD6	PD6 + PD5	EE3-Sunset, EE2-Violet, EE5-Sky	C3-Dark
A3	C3	526 (5M3)	E6-Dark	2x EO3 + 1x EE4 Fog	526 (5M3)	OD6 + PD5	PD5 + PD6	EE3-Sunset, EE2-Violet, EE5-Sky	C3-Dark

The product shades referred to in our comparative list are Dentsply Sirona's closest approximation to the referenced VITA 3D-Master shades. For a few highly chromatic shades, it is recommended to mix OD and PD ceramic masses or either of them with a DE ceramic mass.

## Framework preparation

### Zirconia

For more detailed information regarding fabrication of Cercon and other compatible frameworks, please consult the Cercon Directions for Use (DFU).

### Celtra Press

For more detailed information regarding fabrication of Celtra and other compatible frameworks, please consult the Celtra Directions for Use (DFU).

**Celtra Press – minimum framework wall thickness/veneer thickness (mm)\***

Technology	Area		Inlays	Onlays/ Table- tops	Veneers	Anterior crowns	Poste- rior crowns	Anterior bridges	Posterior bridges
								Connector cross-section 16 mm <sup>2</sup>	
Staining technique	Framework wall thickness (fully contoured)	Full-arch	1.0 ≥ isthmus width	1.5	0.6	1.2	1.5	1.2	1.5
		incisal/occlusal	1.5	1.5	0.6	1.5	1.5	1.5	1.5
Cut-back	Framework wall thickness	Full-arch	–	–	0.6	1.2	1.5	1.2	1.5
		incisal/occlusal	–	–	0.4	0.8	0.8	0.8	0.8
	Veneering (thickness)	–	–	0.4	0.7	0.7	0.7	0.7	
Layering technique	Framework wall thickness	full-arch/incisal/ occlusal	–	–	–	0.8	0.8	0.8	0.8
	Veneering (thickness)	–	–	–	0.4-0.7	0.7	0.7	0.7	

\*Connector for 3-unit bridge should have a minimum of a cross-sectional area of 16 mm<sup>2</sup>

\*\* For other Lithium Disilicate frameworks use the appropriate recommendations for minimum wall thickness/veneer thickness (mm), as per the applicable DFUs.

### PowerFire (Celtra Press frameworks only)

Use 50-micron sized aluminum oxide at 20 psi pressure and lightly blast the exterior surface of the ceramic restoration. Care must be taken not to harm the margins.

Use a steam cleaner to clean the surfaces or put the restoration in distilled water and place in an ultrasonic cleaner for 10 minutes.

PowerFire is a firing program that is carried out before the first ceramic firing of the veneering porcelain. PowerFire increases the flexural strength of the Celtra Press restoration to >500 MPa. After PowerFire, **1** blasting must be avoided as it will reduce the strength of the restoration.

#### Helpful Tips:

- To avoid tensions in the Celtra Ceram veneering ceramic, it is necessary to round off any sharp corners and edges of the frame.
- Irrigation will generally not be required when working on the material with diamond cutters.
- Water cooling is recommended when using high-speed rotary instruments.
- To avoid overheating of the framework material, do not use a high grinding pressure.

**Composite dies**

The **Dentsply Sirona Die Material** is designed to mimic the actual shade of the patient's prepared tooth. When this material is placed inside the pressed Celtra crown, it will assist in accurate shade reproductions.

The dentist should take a prep shade of the tooth being restored for the laboratories reference with the Dentsply Sirona Prep Guide.

If the dentist did not take the shade of the prepared tooth, the composite die material guide below may be used to verify the final shade. Select the proper composite die material from the appropriate table.

- 1) Apply the Dentsply Sirona Prosthetics Die Release to the inside of the ceramic restoration and allow it to dry.
- 2) Place Dentsply Sirona Die Release in the inside of the restoration, then, a small amount of the Dentsply Sirona composite die material on the inside of the restoration. Pack the material to remove any voids. Immediately push a dowel pin into the uncured composite die material. Remove any excess composite from the margin area.

- 3) Light cure the composite for 1–2 minutes using a hand held light curing unit or the Triad 2000 curing unit from Dentsply Sirona.
- 4) Remove the composite die material from the restoration and carefully clean it using a steam cleaner or in distilled water in an ultrasonic cleaner for 10 minutes.

Due to the high translucency of Celtra Press, the influence of the die shade on the shade of the restoration must be taken into account. The aesthetic result is also influenced by the color of the adhesive material. Using the supplied light-curing die material, the dental technician has the ability to map the shade information supplied by the dentist to a control die to replicate information about the oral situation in the shade reproduction. The aim is to simulate the shade of the prepared tooth (follow the working instructions).

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**Celtra® Ceram processing**

**Note:** When firing a Celtra® Press All-Ceramic restoration it is important to **use only ceramic/porcelain** type pins/pegs or to place the object directly on the firing pad, to prevent issues with Celtra® Press during porcelain and glaze firings. When other than recommended type of pins/pegs are used, while initial results with some firing pins may appear acceptable, internal stress can compromise long term success. Do not fill the entire restoration with refractory putty. Doing so could cause cracking of restoration.

### Framework-porcelain ratio standards\* for Celtra Press restorations.

	Veneer			Crown & Bridges including 2 <sup>nd</sup> premolar						
	0.8	1.0	1.1	1.2	1.5	1.7	2.0	2.2	2.4	2.8
Overall thickness of restoration (mm)										
Minimum framework thickness (mm)	<b>0.4</b>	<b>0.5</b>	0.6	0.8	0.8	0.9	1.1	1.2	1.3	1.5
Maximum layer porcelain thickness (mm)	<b>0.4</b>	<b>0.5</b>	0.5	0.4	0.7	0.8	0.9	1.0	1.1	1.3

\*The strength of the veneering must not exceed 2.0 mm at any point.

- The minimum framework wall thickness will always be based on the total thickness of the restoration.
- The thickness ratio of the framework wall to the ceramic layer must be at least 1 : 1 to ensure framework stability and esthetics.

#### Cut-back Technique

The cut-back framework is complemented in the incisal or occlusal areas using incisal materials. This produces high-quality esthetic restorations in only a few steps. The exact tooth shade is ensured after applying the incisal material.

Additional dentin materials need not be applied to Celtra Press restorations.

#### Layering Technique

By applying the individual layering technique, one can use the broad range of dentin and enamel powders **2**, **3**, **4**, **10** & **11** giving you the ability to create natural and highly sophisticated restorations.

If the structure requires dentin proceed with the build-up, aesthetic effects can be achieved in the enamel area using enamel opal **5**, **11** & **15** and enamel effect powders **3**, **4**, **10** & **11**. For cervical, occlusal fossa and body areas use EE1 Sunrise **7** & **13** or EE3 Sunset to enhance chromatic effects. Opal effect powders can be applied as an overlay or as needed **5**, **11** & **15** to further enhance individualization effects.

Three unit bridge after first firing **16**.

#### Staining / Overglazing

The Dentsply Sirona Universal Stain and Glaze system (available separately) is recommended for finishing full contour and veneered restorations. This system is especially matched to the Celtra Press framework and Celtra Ceram veneering system.

Place a small amount of the enamel stain or glaze on the palette. If necessary, mix the stain and glaze with the Dentsply Sirona Stain and Glaze Liquid to achieve a creamy consistency and apply the mixture to the porcelain surface.

Additional Universal Stains may be applied for individual characterizations **17** & **18**.

Remove the die material from the restoration. To ensure proper fit, remove excess glaze from the interior of the crown as well as the internal margin areas.

Fire the crown according to recommended firing cycles **19**.

**Note:** For best results, be sure to mix thoroughly the stain and overglaze prior to use. The pigment and liquid can sometimes become separated inside the jar over time.

**Note:** If a higher sheen is desirable, either raise the high firing temperature 10°C or use an additional 30 second hold time at the high temperature.

## General firing recommendations – Cut-back and Layering Technique

### Celtra Press Framework

	Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
<b>PowerFire:</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
<b>Framework only</b>	0:00	1:00	400	1:00	55	760	0	0	0:00	2:00	0:00
<b>First firing:</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
<b>Dentin &amp; Enamel</b>	2:00	2:00	400	2:00	55	770	400	770	1:00	1:00	5:00
<b>Second firing:</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
<b>Dentin &amp; Enamel</b>	2:00	2:00	400	2:00	55	760	400	760	1:00	1:00	5:00
<b>Glaze Firing</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
	2:00	2:00	400	2:00	55	750	0	0	0	2:00	5:00
<b>Add-on (with and after glaze firing)</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
	2:00	2:00	400	2:00	55	750	400	750	1:00	1:00	5:00

### Zirconia Framework

	Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
<b>First firing:</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
<b>Dentin &amp; Enamel</b>	2:00	2:00	400	2:00	55	780	400	780	1:00	1:00	0:00
<b>Second firing:</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
<b>Dentin &amp; Enamel</b>	2:00	2:00	400	2:00	55	770	400	770	1:00	1:00	0:00
<b>Glaze Firing</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
	2:00	2:00	400	2:00	55	750	0	0	0:00	2:00	6:00
<b>Add-on (with and after glaze firing)</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
	2:00	2:00	400	2:00	55	750	400	750	1:00	1:00	6:00

## Full Contour Technique

### Celtra Press Framework

	Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
<b>Power firing incl. Glaze (Celtra Press framework only) - 1st firing</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
	2:00	2:00	400	2:00	55	760	0	0	0	2:00	5:00
<b>Glaze – 2nd firing</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
	2:00	2:00	400	2:00	55	750	0	0	0	2:00	5:00
<b>Add-on with 1st glaze firing</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
	2:00	2:00	400	2:00	55	760	400	760	1:00	1:00	5:00
<b>Add-on after glaze firing</b>	min	min	°C	min	°C/min	°C	°C	°C	min	min	min
	2:00	2:00	400	2:00	55	750	400	750	1:00	1:00	5:00

## Full Contour Technique

### Zirconia Framework

	Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
<b>Glaze Firing</b>	min 2:00	min 2:00	°C 400	min 2:00	°C/min 55	°C 760	°C 0	°C 0	min 0:00	min 2:00	min 6:00
<b>Add-on (with and after glaze firing)</b>	min 2:00	min 2:00	°C 400	min 2:00	°C/min 55	°C 760	°C 400	°C 760	min 1:00	min 1:00	min 6:00

- Note:**
1. Slow cooling is mandatory; this includes correction firings of restorations after try-in.
  2. Firing temperatures must be adapted to the number of units fired in the same cycle.
    - a. 5 to 9 units require an increase by 5 °C to 10 °C;
    - b. 10 or more units require an increase by 10 °C to 20 °C.

The values indicated here are recommended values and serve only as a guidance. Deviations of the firing results are possible. The firing results depend on the respective output per oven and are due to the manufacturer and to age. Therefore the recommended values must be individually adapted at each firing.

We recommend a test fire to control the oven. All indications have been carefully elaborated and tested by us, but they are passed on without any guarantee.

For up-to-date firing recommendations please visit [celtra-dentsplysirona.com](http://celtra-dentsplysirona.com).

## Cementation

### Preparation of the Celtra restoration

- Clean the restoration with a steam cleaner, in ultrasonic bath or with alcohol.
- Apply 5% – 9% hydrofluoric acid etching gel (Available separately, see manufacturer's complete Directions for Use) to the interior of the restoration only and allow to soak for 30 seconds.
- CAUTION: Follow manufacturer's precautions. Do not allow tissue or eyes to come into contact with the acid!
- Remove the hydrofluoric acid as per the manufacturer's instructions.
- Dry the restoration in an air stream. It is recommended to silanize the etched surfaces immediately.
- At chairside, apply silane only to those surfaces required for adhesive cementing.
- Allow to soak for 60 seconds. If the silane layer is no longer liquid, add more silane. Blow-dry in a powerful air stream. (Recommended material: Calibra® Silane Coupling Agent, available separately, see complete Directions for Use).

### Cementing

Depending on the indication for Celtra® Press restorations a self-adhesive or fully-adhesive cementation can be chosen. Compatible time-proven adhesive cementing materials are available as part of the Dentsply Sirona range of products. Alternatively, crowns and bridges can also be fixed with glass ionomer cement. Cements are available separately.

	Self-adhesive	Fully adhesive	Glass-ionomer
Inlays	R	HR	–
Onlays	R	HR	–
Veneers	–	HR	–
Crowns	HR	HR	R
Bridge	R	HR	R

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