Etchant Gel S

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Note: Gloves, safety goggles (for both doctor and patient), and mouth guards are important for safety reasons! Avoid bringing Etchant Gel into contact with gingiva or adjacent teeth. The needle tip is single-use only. To prevent cross-

contamination discard it after use.

1. Cleaning teeth

Clean the tooth which is to be treated using a fluoride-free cleaning paste.

| **2. Keeping work dry (rubber dam)** Optimal results can only be achieved if the work area is dry. Use of a dental dam is recommended. |
| **3. Preparation of the Cavity** Cause acid burns. Done in accordance with the principles of adhesive restoration methods. For pulp protection, a small amount of hard setting calcium hydroxide cements should be placed close to the pulp. |
| **4. Application** Apply Etchant Gel S directly from the applicator needle on the surface to be treated. Allow etchant to remain in place for 30–60 s on enamel (depending on preparation) and for 15 s on dentin. Afterwards, rinse thoroughly with water (approx. 20 s), until all traces of Etchant Gel have been removed. Gently dry with pressurized oil-free air, taking care not to overdry the dentine. Proceed as indicated in bonding manufacturer’s instructions. Etched enamel has a chalky white appearance. In case of doubt repeat etching procedure. |
| ***Note: Use fluoride on enamel which has been etched and is not covered by composite.*** |

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