Zeta 7 Solution

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EN - Concentrate disinfectant liquid for impressions

Zeta 7 Solution is an aldehyde-free concentrated disinfectant with a wide range of action prepared and tested in compliance with new harmonised European disinfection regulations. Zeta 7 Solution has been designed for high-level, rapid disinfection of impressions (addition and condensation), in silicone, alginate, polyether and polysulphur: the impressions remain unaltered, have no significant dimensional change and are compatible with plaster.

COMPOSITION: 100 g of Zeta 7 Solution contains 7.7 g dimethyl-didecyl-ammonium chloride, 15 g phenoxyethanol, non-ionic surfactants, additives, auxiliaries and water to 100 g.

USE: Immediately after taking the impression, rise the impression under running water for 30 seconds. Inspect the impression to make sure that all visible contamination has been removed. In case of visible contamination repeat rinsing. Place the impression for 10 minutes in a 1% solution bath of Zeta 7. Use the dispenser cap to measure 10 ml of Zeta 7 Solution to be added to every litre of water. Remove the impression from the bath and rinse it thoroughly. Follow recommended immersion times as suggested. Replace the disinfectant solution every day, and in any case never disinfect more than 40 impressions with the same solution. In the day, replace solutions immediately if they are visible contaminated (blood, saliva or organic tissue) and/or if the colour/odour of the disinfectant solution changes, as this may interfere with its disinfectant action. Unused disinfectant solutions will remain stable for 7 days if stored in a sealed plastic container (HDPE, GPPS, ABS, PP and PMMA) away from light and heat.

Efficacy (microorganisms)	Reference standard	Concentration	Reaction time (minutes)
Bactericidal (S. aureus, P. aeruginosa, E. hirae)	EN 13727	1%	10
Yeasticidal (C. albicans)	EN 13624	1%	10
Tuberculocidal (M. terrae)	EN 14348	1%	10
Virucidal against enveloped viruses including HIV, HBV, HCV	EN 14476	1%	10

WARNINGS, PRECAUTIONS AND CONSERVATION CONDITIONS: this product has been designed for use by qualified operators with knowledge of the required safety standards. Use protective gloves. **Operate in well-ventilated areas.** Store the closed container in a cool dry place. Prevent Zeta 7 Solution from coming into contact with other disinfectants or cleaners. Any slight difference in the perfume is due to the presence of natural fragrances and does not jeopardize quality of the product.

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