

Ultimate Etchant Gel



coltène®
ETCHANT GEL S

coltène®
ETCHANT 15

Efficient and Safe

coltène® ETCHANT GEL S

coltène® ETCHANT 15

Precise Placement

The thin, bendable tip allows you to place Coltène® ETCHANT GEL exactly where you need it, even in hard to reach areas.



Smooth Consistency

Flows easily into place without dripping or slumping.

Easy Removal

Rinses off easily. ETCHANT GEL's bright colors provide clear indication of thorough rinsing.

Ultimate ETCHANT GEL

ETCHANT GEL does not dry out, maintaining its smooth consistency over months of repeated use.



Ordering information

Coltène® ETCHANT GEL is available in 35 % (blue) or 15 % (green) phosphoric acid gel; in syringes or bottles.

7373 ETCHANT GEL S Kit	35 %
Syringes	2 × 2.5 ml
Application Needles	20 pcs.
7886 ETCHANT GEL S Refill	35 %
Syringes	4 × 1.2 ml
Application Needles	20 pcs.
7374 ETCHANT GEL	35 %
Bottle	1 × 5 ml
7878 ETCHANT 15 Kit	15 %
Syringes	4 × 1.2 ml
Application Needles	20 pcs.
7876 Application Needles	black
	20 pcs.

www.coltenewhaledent.com
© 2007 Coltène/Whaledent AG

Coltène/Whaledent Ltd.
The President Suite-A
Kendal House, Burgess Hill
Victoria Way
West Sussex, RH15 9NF/U.K.
Tel. +44 (0)1444 235486
Fax +44 (0)1444 870640
hq@coltenewhaledent.co.uk

Coltène/Whaledent AG
Feldwiesenstrasse 20
9450 Altstätten/Switzerland
Tel. +41 (0)71 757 53 00
Fax +41 (0)71 757 53 01
office@coltenewhaledent.ch
www.coltenewhaledent.com

Coltène/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, Ohio 44223/USA
Tel. +1 330 916 8800
Fax +1 330 916 7077
productinfo@coltenewhaledent.com
www.coltenewhaledent.com

coltène
whaledent