

I AM  
VERSATILE

AIR N GO  
EASY



**All clinical applications  
are easy with the  
all-in-one air polisher**

**ACTEON**

# SIMPLE AIR POLISHING

## AIR N GO

EASY

### DUAL FUNCTIONS

- "SUPRA" function for supragingival polishing and stain removal
- "PERIO" function for subgingival polishing, periodontal maintenance and peri-implant treatments



**A single handpiece to cover a range of care and treatment**

**AIR-N-GO® easy** combines essential functions and reliability:

- Fast start-up: direct connection to the high speed handpiece connector
- 360° swivel handpiece for fast and accurate treatment
- Parts can be dismantled for easy cleaning

### 4 Nozzles for Treatment

- Sodium bicarbonate "CLASSIC" powders
- Glycine "PERIO" powder



# A COMPLETE SOLUTION FOR IMPLANT MAINTENANCE

NEWTRON<sup>®</sup>



NEWTRON



implantProtect  
Pure Titanium



AIR N GO  
EASY



IMPLANT PROTECT TIP KIT

AIR-N-GO<sup>®</sup> + PERIO  
EASY + EASY KIT

DETECT

DIAGNOSE

TREAT

Contact ACTEON for a Detailed Implant Maintenance Protocol

# PREVENT & FIGHT DISEASE

**AIR-N-GO® nozzles provide continuous, accurate and controlled spray. These nozzles offer power and efficiency for the removal of supra and subgingival biofilm.**

## AIR-N-GO® SUPRA NOZZLE (standard)



*Clinical case Dr. P. ROUAS (France)*

- Suitable for any supra-gingival prophylactic treatment in combination with **AIR-N-GO® "CLASSIC"** powder


## AIR-N-GO® PERIO easy NOZZLE (optional Perio easy kit)



*Clinical case Dr. C. LIM (U.K.)*

- Innovative design in a "shoehorn" shape providing easy sub-gingival insertion
- For non-surgical, sub-gingival treatment of teeth and implants

## AIR-N-GO® PERIO MAINTENANCE NOZZLE (optional)



*Clinical case Dr. C. LIM (U.K.)*

- Reduced consumption of **AIR-N-GO® "PERIO"** fine powder particles supragingivally
- For use on implants and appliances

## AIR-N-GO® PERIO NOZZLE (optional)



*Clinical case Dr. A. PATEL (U.K.)*

- Sub-gingival use for surgical and non-surgical treatments (tunnelling and open flaps)
- Powerful double lateral spray
- For use on implants



# CUSTOMIZED SOLUTIONS FOR YOUR CLINICAL NEEDS

## AIR-N-GO® "CLASSIC" POWDERS

(Sodium bicarbonate)  
(≈ 76µm average grain size)


### Complete range of active supragingival powders

- Cleaning provides instant results
- Remove stains and bacterial plaque, even in areas that are difficult to access such as pits and fissures
- Better patient comfort: gentle and efficient with limited bleeding during treatment



**5 fresh flavors based on 100% natural aromas or essential oils:  
Neutral, Cola, Raspberry, Peppermint  
and Lemon**

## AIR-N-GO® "PERIO" POWDER



AIR-N-GO® "PERIO" (x1000)  
Source: ICMCB CNRS (Institute of Condensed Matter Chemistry - French National Centre for Scientific Research) S.E.M. photos



**SODIUM  
FREE!**

*Natural glycine-based  
(≈ 25µm average grain size)*

- Low abrasiveness preserving the treated surfaces and the most fragile anatomical areas
- Clinically proven to decrease biofilm subgingivally
- Safe to use on any implant surface

# AIR N GO

EASY



GO SUPRA

Teeth and Implants

GO PERIO

Nozzles	Powders	Indications
 <p>AIR-N-GO® easy SUPRA 120° Supplied as standard</p>	 <p>"CLASSIC" Sodium bicarbonate-based Average particle size: 76µm FLAVORS: Neutral, Cola, Raspberry, Peppermint and Lemon</p>	<ul style="list-style-type: none"> <li>• Finishing after scaling</li> <li>• Preparation for whitening</li> <li>• Cleaning of interproximal areas, fissures and crevices</li> <li>• Preparation of the tooth surface before etching</li> <li>• Cleaning of orthodontic brackets</li> <li>• Removal of supragingival biofilm, plaque, stains and temporary cement residues</li> </ul>
 <p>AIR-N-GO® PERIO easy Supplied with the PERIO Kit (optional)</p>	 <p><b>SODIUM FREE!</b></p>	<ul style="list-style-type: none"> <li>• Removal of the sub-gingival biofilm</li> <li>• Treatment of periodontal disease and peri-implantitis</li> </ul>
 <p>AIR-N-GO® PERIO (optional)</p>		<p>Treatment of periodontal disease and peri-implantitis</p>
 <p>AIR-N-GO® PERIO MAINTENANCE (optional)</p>	<p>"PERIO" Glycine-based Average particle size: 25µm FLAVOR : Neutral</p>	<p>Periodontal maintenance</p>

F0337/ May 2019