

A fencer in white gear is shown in a dynamic pose, lunging forward with a foil. The background is a light, neutral color.

MAILLEFER

ProTaper Next[®]

Flexible performance



WE
KNOW
ENDO.



Dentsply
Sirona



GLIDE PATH



PG

ProGlider®

SHAPING



X1

X2

X3

X4

X5

ProTaper Next®

Motor Settings:
300 rpm / 2 - 5.2 Ncm

ProGlider[®] Glide Path

PG



016

2%*

ProTaper Next[®] Shaping

X1



017

4%*

X2



025

6%*

X3



030

7%*

X4



040

6%*

X5



050

6%*

*variable taper

Choose your obturation solutions



ProTaper Next®
paper points



Protaper Next®
Conform Fit™ gutta-percha



GuttaCore® for
ProTaper Next®

INL/ B EN PTNO TIP 000 / 04 / 2017 -
updated 01/2021



Mailefer Instruments Holding Sàrl
Chemin du Verger 3 - 1338 Ballaigues - Switzerland
Phone: +41 (0)21 843 92 92 - endo@dentsplysirona.com
www.dentsplysirona.com

