

Designed to Provide a
Better Solution for
Crown and Bridge Preparation

GREAT WHITE[®]
ULTRA
CROWN & BRIDGE PREPARATION CARBIDE BURS





simply better dentistry™
WHITE[®]

Faster Bulk Reductions with a Single Instrument

Cut up to 40% faster with significantly less effort than diamonds*

Great White® Ultra provides high-speed performance and cuts restorative materials faster, enabling transition between different substrates with minimal clogging and no drag.

This results in less operative pressure and easier instrument advancement.



The Great White® Ultra:

- Cuts quickly and smoothly through enamel
- Negotiates amalgam and other restorative materials with minimal clogging and no “drag or stalling”
- Allows completion of the bulk reduction task in the crown preparation with a single instrument
- More efficient –less chair time
- Helps provide precise margins (Please reference SEM photograph)

GREAT WHITE[®] ULTRA

CROWN & BRIDGE PREPARATION CARBIDE BURS

FOOTBALL



GWU 379-023

5 PK
HEAD LENGTH (mm) 4.2
MAX HEAD DIAMETER (mm) 2.3

SHORT FLAT END TAPER



GWU 845KR-016

18195-5
4.0
1.6



GWU 845KR-018

18196-5
4.0
1.8



GWU 845KR-025

18197-5
4.0
2.5

FLAT END TAPER



GWU 847-016

5 PK
HEAD LENGTH (mm) 8.0
MAX HEAD DIAMETER (mm) 1.6



GWU 847-018

18191-5
8.0
1.8



GWU 847-020

18192-5
8.0
2.0

FLAT END TAPER SHORT SHANK



GWU 847-016S

5 PK
HEAD LENGTH (mm) 8.0
MAX HEAD DIAMETER (mm) 1.6



GWU 847-018S

18187-5
8.0
1.8



GWU 847-020S

18185-5
8.0
2.0

ROUND END TAPER



GWU 855-025

5 PK
HEAD LENGTH (mm) 7.0
MAX HEAD DIAMETER (mm) 2.5



GWU 856-016

18198-5
8.0
1.6



GWU 856-018

18188-5
8.0
1.8



GWU 856-020

18189-5
8.0
2.0



GWU856-023

18190-5
8.0
2.3

ROUND END TAPER SHORT SHANK



GWU 856-016S

5 PK
HEAD LENGTH (mm) 8.0
MAX HEAD DIAMETER (mm) 1.6



GWU 856-018S


18183-5
8.0
1.8



GWU 856-020S

18184-5
8.0
2.0

Crown and Bridge Preparation Kit



Images
Not Shown
to Scale

All-In-One Kit

Designed for the efficient atraumatic removal of failed amalgams, composites, zirconia and PFM Crowns.

Kit Reorder No. 18150

Kit Contents:

2-Great White[®] Ultra #856-016
2-Great White[®] Ultra #856-018
2-Great White[®] Ultra #856-020
1-Great White[®] Ultra #379-023
2-Great White[®] Ultra #847-016
2-Great White[®] Ultra #847-018
1-Autoclavable Bur Block

The Great White[®] Ultra:

- Cuts quickly and smoothly through enamel
- Negotiates amalgam and other restorative materials with minimal clogging and no "drag or stalling"
- Allows completion of the bulk reduction task in the crown preparation with a single instrument
- More efficient—less chair time
- Helps provide precise margins (Please reference SEM photograph)

"The Great White[®] Ultra is a versatile bur designed for crown and bridge. It cuts smoothly and efficiently through dentin, enamel, amalgam, composite and cast metal restorations. The unique blade geometry ensures smooth cutting action even when moving from dentin to amalgam or cast metal restorations."

— Steven Abrams, DDS

Great White® Ultra Brochure - Rev.2



6 4 6 7 4

800.535.2877 | www.sswhitedental.com
1145 Towbin Avenue Lakewood, New Jersey 08701