



PRODUCT CATALOGUE

Amber Rosetta Zirtooth



CAD/CAM blocks

Compatible block product line-up based on the material loading, machinability, and wet/dry status of each user's CAD/CAM equipment.



	Rosetta [®] SM		Amber [®] Mill	Amber [®] Mill DIRECT	Amber [®] Mill 
	Blocks	Disk			
	HASS's 1st generation	Ø98 Rosetta SM Disk	Freedom of transparency according to heat treatment. HT → MT → LT → MO 	No need to heat treatment Functional gradient block	Hybrid ceramic block (glass + resin)
					
Indications	Inlays				
	Onlays				
	Veneers				
	Anterior Single Crowns				
	Premolar Single Crowns				
Posterior Single Crowns				Possible	Possible
	3-Unit Bridge (up to the 2nd Premolar)			Impossible	Impossible
Flexural Strength (MPa)	440 MPa	440 MPa	450 ± 42 MPa	avg. 355 MPa 300±15 MPa – 380±25 MPa Incisal → Cervical	200 MPa
Do need firing					Don't need firing
Translucency & Shade	HT	W1 W2 W3 W4 A1 A2 A3 A3.5 B1 B2 B3 B4 C1 C2 C3 C4 D2 D3 D4	W1 W2 W3 W4 A1 A2 A3 A3.5 B1 B2 B3 B4 C1 C2 C3 C4 D2 D3 D4	W1 W2 A1 A2 A3 A3.5 B1 B2 C2 D2	A1 A2 A3 B1
	MT				
	LT	W1 W2 W3 W4 A1 A2 A3 A3.5 B1 B2 B3 B4 C1 C2 C3 C4 D2 D3 D4	W1 W2 W3 W4 A1 A2 A3 A3.5 B1 B2 B3 B4 C1 C2 C3 C4 D2 D3 D4		W2 A1 A2 A3 A3.5 B1
	MO	MO0 MO1 MO2 MO3 MO4	MO0 MO1 MO2 MO3 MO4		

W1 W2 W3 W4
A1 A2 A3 A3.5
B1 B2 B3 B4
C1 C2 C3 C4
D2 D3 D4

The transparency
can be adjusted
according to the
heat treatment,
so only the Shade
is selected.