


Scanning without powder

Registrado® Scan



Registrado® Scan

SCANNABLE ADDITION-CURING SILICONE FOR BITE
REGISTRATION, HARD, FAST SETTING

Registrado® Scan

BITE REGISTRATION MATERIAL, SCANNABLE WITHOUT POWDERING – A-SILICONE, HARD, FAST SETTING

With Registrado Scan VOCO is leading its proven product line Registrado into the digital age. Collecting 3D data when recording an occlusion requires a material suitable for digital workflow to ensure that the occlusion is reproduced with precision.

Thanks to its special visual properties, optimised in brightness and contrast, Registrado Scan satisfies the strict requirements of digitalisation. Its visual properties allow trouble-free scanning of the occlusal record without using scanning powder. The occlusal record can be digitalised very quickly after its removal from the mouth without any additional steps.



Digitisable occlusion record with special visual properties

Registrado Scan is based on the proven Registrado quality. The material covers a wide range of indications, from practical working time and consistency appropriate for the indications to good final hardness, perfect reproducibility and dimensional stability.

You choose: With Registrado Scan you can work conventionally as usual. After recording an occlusion you decide if you are going to continue your work digitally.

Indications

- Occlusal bite registration
- Fixing of face bow registrations
- Location material for intraoral pin tracing templates

Advantages

- Scanning without powder
- Precise reproduction of occlusion
- High Shore A hardness [90]
- Short intraoral setting time [45 sec.]
- Good cutting properties
- Analog or digital

Presentation

- REF 2038 Cartridge 50 ml, accessories
- REF 2158 Dispenser – type 2
- REF 2196 Mixing tips type 16, 50 pcs.
- REF 2140 Intraoral tips type 2, 50 pcs.



VOCO GmbH
Anton-Flettner-Straße 1-3
27472 Cuxhaven
Germany

Freecall: 00 800 44 444 555
Fax: +49 (0) 4721-719-140

info@voco.com
www.voco.com

Available from: