



# KENDA NOBILIS

1-STEP DIAMOND POLISHING SYSTEM

- Polish porcelain and composites with the same polisher
- For a silky luster in only 1 step
- Economical diamond polisher ideal for posterior restorations
- Flexible, adapting to tooth shape
- Recommendation: For anterior restorations, use KENDA Maximus or KENDA Unicus polishers subsequently for ultimate high shine polishing



SET REF 0300.006A



<b>REF</b> reorder	0301	0302	0303	0305	0306	0308	0309	0310	0311
<b>ISO</b> 1: 858 204...	243 513 055	254 513 050	304 513 100	030 513 060	243 513 030	243 513 050	303 513 100	030 513 090	030 513 060
$\varnothing \times L$ mm	5,4 x 13	5 x 8	10 x 2	6,5 x 10	3 x 7	5 x 10	10 x 1	9 x 8,5	6 x 7,5
$\mu$ grit	fine		$\mu$ grit		color		rec.	max.	
	fine		violet		violet		7500 - 10 000 RPM	10 000 RPM	
Refill Box of:	N Pressure:		Autoclavable:		Shank:		Medical device:		
6, 12, 25, 100	0,3 - 2N		135°C, 275°F		RA, stainless steel		Class IIa		

