



# GC New Metal Strips

Metal Finishing Strips

## Epitex

Synthetic/Plastic Finishing and Polishing strips




### New Metal Strips from GC

Color-coded metal strips, electro-plated with aluminium-oxide particles for approximal contouring, finishing and polishing of all restorative filling materials.

### Epitex from GC

Strips for contouring, finishing and polishing in interproximal spaces, suitable for glass ionomers and composites.

**GC**



## GC New Metal Strips

- Advantages**
- › easy to handle, flexible and long-lasting
  - › ultra-thin, yet very strong
  - › can be sterilized or autoclaved repeatedly
  - › various grades and widths (2,6-3,3-4,0 mm)

- Indications**
- › N° 200 (red) - preparation
  - › N° 300 (blue) - final preparation
  - › N° 600 (green) - finishing
  - › N° 1.000 (yellow) - polishing



## GC EpiteX

- Advantages**
- › Extremely thin and flexible material
  - › Excellent adaptability even between tight interproximal contacts
  - › Very tear resistant
  - › Abrasive particles are firmly embedded into the strip surface
  - › Four grades
  - › Handy strip dispenser
  - › Also available as matrix strips

- Indications**
- › N° 200 (red) - preparation
  - › N° 300 (blue) - final preparation
  - › N° 600 (green) - finishing
  - › N° 1.000 (yellow) - polishing



Assortment pack

GC New Metal Strips	
<b>Assortment pack</b>	
000267	GC New Metal Strips assortment Content 12 pcs. (1 pcs. of each width and grain)
<b>Content: 12 pcs. (4pcs. each of widths 2.6mm, 3.3mm, 4.0mm)</b>	
000251	GC New Metal Strips Red Mesh N° 200 (140µm Rough Preparation)
000255	GC New Metal Strips Blue Mesh N° 300 (90µm Preparation)
000259	GC New Metal Strips Green Mesh N° 600 (70µm Finishing)
000263	GC New Metal Strips Yellow Mesh N° 1000 (50µm Polishing)
<b>6mm Refill packs: (12pcs. of 6mm wide)</b>	
000268	GC New Metal Strips Red Mesh N° 200 (140µm Rough Preparation)
000269	GC New Metal Strips Blue Mesh N° 300 (90µm Preparation)
000270	GC New Metal Strips Green Mesh N° 600 (70µm Finishing)
000271	GC New Metal Strips Yellow Mesh N° 1000 (50µm Polishing)
GC EpiteX	
<b>Assortment pack</b>	
000404	GC EpiteX, Starter Kit
<b>Content: Stand with dispenser containing all available grains, 10m each (coarse-blue, medium-green, fine-grey, extra fine-red)</b>	
<b>Unit refill packs (Content: 10m)</b>	
000409	GC EpiteX Translucent matrix, 10 m
000405	GC EpiteX Coarse (blue), 10 m
000407	GC EpiteX Medium (green) 10 m
000406	GC EpiteX Fine (grey), 10 m
000408	GC EpiteX Extra fine (red), 10 m

### GC EUROPE N.V.

Head Office  
 Researchpark Haasrode-Leuven 1240  
 Interleuvenlaan 33 - B-3001 Leuven  
 Tel. +32.16.74.10.00 - Fax. +32.16.40.48.32  
 info.gce@gc.dental  
<http://www.gceurope.com>

### GC UNITED KINGDOM Ltd.

Coopers Court  
 Newport Pagnell  
 UK-Bucks. MK16 8JS  
 Tel. +44.1908.218.999  
 Fax. +44.1908.218.900  
 info.uk@gc.dental  
<http://uk.gceurope.com>

